



GSA prices effective October 7, 2024  
(based on July 1, 2024, commercial price list)  
Published March 2025

## Laboratory Furniture Co/Struc® System

### General Services Administration Federal Supply Service Authorized Federal Supply Schedule Price List

Online access to contract ordering information, terms and conditions, up-to-date pricing, and the option to create an electronic delivery order is available through GSA Advantage!®, a menu-driven database system. The Internet address for GSA Advantage! is [www.gsaadvantage.gov](http://www.gsaadvantage.gov).

Office Furniture  
FSC/PSC 7110/7125  
Contract Number: GS-03F-036DA  
UEI: Q2K3MSZ843D8  
Contract Period: December 21, 2015—December 20, 2025

MillerKnoll, Inc.  
855 East Main Avenue  
Zeeland, Michigan 49464  
(616) 654 3000 Phone  
(616) 654 8278 Fax  
[www.HermanMiller.com/government](http://www.HermanMiller.com/government)

Business Size: Large



# Contract Information

## 1a. Special Item Numbers (SIN)

SIN 33721

Furniture Systems, Computer Furniture, Filing and Storage, Tables and Accessories, Upholstered Seating and Multi-Purpose Seating, Project Management, Reconfiguration and Relocation Services, Design/Layout and Installation Services

SIN 33721P

Packaged Offices

SIN 337127

Modular Laboratory Furniture Systems

SIN 339113H

Healthcare Furniture

SIN OLM

Order-Level Materials (OLMs)

This SIN is only for products and/or services used in direct support of the purchase of new furniture. It requires approval from the ordering Contract Officer. Please contact your local MillerKnoll Inc. representative to discuss the appropriate use of this SIN.

SIN 532289

Leasing

Contact MillerKnoll Inc participating government dealer for current leasing information.

## 1b. Lowest Priced Model

SIN 33721:	BVP3.S	\$7.27
SIN 33127:	TW113	\$5.81
SIN 339113H:	CR900NR	\$483.80

## 1c. Service Rates

Project Management

Furniture Project Manager: \$64.75 per hour  
Furniture Project Coordinator: \$49.00 per hour  
Build-Out Project Manager: \$80.00 per hour

Reconfiguration and Relocation Services

Not to exceed \$55 per hour.

See \* Additional Notes for Installation & Reconfiguration and Relocation Services.

Design/Layout

Not to exceed \$65 per hour.

Installation Services

Not to exceed \$55 per hour (\$82.50 per hour for after-hours Installation).

See \* Additional Notes for Installation & Reconfiguration and Relocation Services.

## 2. Maximum Order

SIN 33721, 337127, NEW, OLM

\$250,000 net product value

SIN 33721P

\$250,000 net value

SIN 339113H

\$500,000 net product value

## 3. Minimum Order

\$100 net

## 4. Geographic Coverage

Continental U.S.

Orders for Alaska, Hawaii, and Puerto Rico will be shipped F.O.B. inland carrier point of exportation.

## 5. Point of Production

Zeeland, Ottawa County, Michigan

Spring Lake, Ottawa County, Michigan

Elland, Yorkshire, United Kingdom

Hildebran, Burke County, North Carolina, USA

## 6. Pricing

Prices shown are list. Discounts include Industrial Funding Fee.

SIN 33721

Standard Lead Time, Single Order Net Purchase \$100 - 250,000, includes approved Herman Miller Options Products.

Product Line	Discount From List
Action Office 1	73.8%
Action Office 2	73.8%
Aeron	64.5%
Ambit Workspace Solutions	74.3%
Asari Chair	53.7%
Aside	64.2%
Bay Work Pod	56.7%
C-Style Overhead Storage	73.8%
Canvas Office Landscape Channel	74.8%
Canvas Office Landscape Dock	74.8%
Canvas Metal Desk	74.8%
Canvas Office Landscape Metal Storage	74.8%
Canvas Office Landscape Private Office	74.8%
Canvas Office Landscape Wood Storage	74.8%
Canvas Office Landscape Vista	74.8%
Canvas Office Landscape Wall	74.8%
Caper	61.5%
Commend Nurses Station	73.3%
Cosm	61.7%
Eames® Aluminum Group	58.7%
Eames Conference & Dining Tables	73.6%
Easton Family	54.91%
Embody	59.7%
Energy Distribution Systems	73.8%
Ethospace System	73.8%
Everywhere™ Tables	73.8%
Exclave	61.7%

Contract Information

continued

Fuld	59.7%
Headway™ Tables	73.8%
Intent Solution	58.7%
Layout Studio®	60.7%
Lighting	73.8%
Lino	63.7%
Meridian Laterals	67.3%
Meridian Pedestals	64.2%
Meridian Towers	62.7%
Meridian Storage	62.7%
Mirra 2	61.7%
Motia Tables	73.8%
NaughtOne	50.0%
Nemschoff	54.91%
Nevi™ Tables	73.8%
Nevi Link	73.8%
OE1 Boundary	72.8%
OE1 Community Elements	61.7%
OE1 Micro Packs	73.8%
OE1 Storage	72.8%
OE1 Tables and Benching	73.8%
Overlay™	62.7%
Plex® Lounge Furniture	59.5%
Pronta Stacking Chair	57.7%
Public Office Landscape	61.7%
Renew™ Tables	73.8%
Renew Link	73.8%
Sayl	62.2%
Setu	59.7%
Swoop Lounge Furniture	60.7%
Thrive Portfolio	62.0%
Thrive Portfolio-CBS Products Extension	62.0%
Tu Lateral Files	74.3%
Tu Pedestals	73.8%
Tu Storage	74.3%
Tu Towers	74.3%
Tu Wood Cases	74.3%
Tu Wood Cubbies	74.3%
Tu Wood Credenzas	74.3%
Tu Wood Pedestals	74.3%
Tu Wood Towers	74.3%
Valor™ Family	54.91%
Verus	63.7%
Zeph	60.2%

Services

Please refer to section 19 for Terms and Conditions of Installation/ Reconfiguration and Relocation Services.

Project Management

Services under this SIN may be used only in conjunction with the purchase of new furniture. Please contact your local MillerKnoll Inc. representative to discuss the available labor categories (listed below). The hourly rates listed are “Not To Exceed” or “NTE” values.

Furniture Project Manager: \$64.75 per hour  
Furniture Project Coordinator: \$49.00 per hour  
Build-Out Project Manager: \$80.00 per hour

Reconfiguration and Relocation Services

Reconfiguration and Relocation Services will be provided by MillerKnoll Inc. or its participating dealer based on a fee not to exceed \$55 per hour.

See \* Additional Notes for Installation & Reconfiguration and Relocation Services.

Design/Layout

Design service will be provided by MillerKnoll Inc. or its participating dealer based on a fee not to exceed \$65 per hour.

Installation Services

Installation Services will be provided by MillerKnoll Inc. or its participating dealer based on a fee not to exceed \$55 per hour (\$82.50 per hour for after-hours Installation).

See \* Additional Notes for Installation & Reconfiguration and Relocation Services.

SIN 337127

Standard Lead Time, Single Order Net Purchase \$100 - 250,000, includes approved Herman Miller Options Products.

Product Line	Discount From List
Action Lab	58.5%
Compass System	58.5%
Co/Struc System	58.5%
Healthcare Carts	18.0%
Mora System	58.5%
Bedside Cabinets	58.5%

SIN 33721P

Packaged Office Furniture - Maximum Order \$250,000 (net).

Includes packaged furniture solutions for customers needing to furnish an office.

SIN 339113H

Standard Lead Time, Single Order Net Purchase \$100 - 500,000, includes approved Herman Miller Options Products.

Product Line	Discount From List
Nemschoff Healthcare Furniture	54.91%

7. Quantity Discounts

Reference item #6.

8. Payment Terms

Net 30 days

9. Foreign Items

None

10. Time of Delivery

- a. 90 day ARO
- b. Expedited delivery: Items listed in the 10-Day section of the catalogs are available for expedited delivery.
- c. Overnight/2-day delivery: Select items may be available for overnight or 2-day delivery. Contact MillerKnoll Inc. or a participating dealer for availability and rates.
- d. Urgent requirements: Agencies may contact contractor or participating dealer to request faster delivery.

# Contract Information

*continued*

## 11. F.O.B.

F.O.B. destination except Alaska, Hawaii, and Puerto Rico, which will be shipped F.O.B. inland carrier point of exportation.

## 12. Ordering Address

a. MillerKnoll Inc.

Government Customer Care 0161  
855 East Main Avenue  
Zeeland MI 49464.

b. For supplies and services, the ordering procedures information on Blanket Purchase Agreements (BPA's) are found in Federal Acquisition Regulation (FAR) 8.405-3.

## 13. Payment Address.

MillerKnoll Inc.  
22764 Network Place  
Chicago, IL 60673-1227

## 14. Warranty

MillerKnoll Inc. commercial warranty applies.

## 15. Export Packing Charge

Prices supplied on request.

## 16. Rental Maintenance and Repair Terms

Not applicable

## 17. Terms and Conditions of Installation/Reconfiguration and Relocation Services

\* Additional Notes for Installation & Reconfiguration and Relocation Services:

FOB Destination (Drop Ship): Price includes product delivery to the site, brought to the tailgate of the truck. The purchaser is responsible for unloading.

Prices effective up to the maximum order by SIN only. Above the MO, pricing shall be negotiated on a case-by-case basis between the client and local dealer.

Standard Services are Subject to the Following Conditions:

Installation will be performed during normal weekday working hours.

Adequate facilities for delivery, unloading, moving and staging/storing the product during the installation process shall be provided.

Service work will not be hindered by other trades.

Electric, heat, and adequate elevator service will be furnished without charge.

The immediate installation area shall be complete and free of debris including the carpet/flooring before installation commences.

Any work requiring a licensed electrician is the responsibility of the buyer.

Additional Charges May Apply For:

Major Metro Markets and any non-ground floor installation.

Major Metro Markets include large population centers and urban environments.

Installation in a clinical/medical environment.

Special restrictions or limits established by local laws, ordinances or the directions of the buyer, including but not limited to restrictions on transportations of materials, street access to the job site and/or dock facilities.

Installations outside of a 50 mile radius of the servicing dealer.

Local Prevailing Wage and/or Union Labor Rates.

Any additional charges shall be quoted by the dealer and approved by the buyer prior to performance of the work.

## 18. Terms and Conditions of Repair Parts

Not applicable

## 19. Service and Distribution Points

For service and distribution points nearest you, call (800) 851 1196.

## 20. Participating Dealers

Access Herman Miller web site for current list of participating dealers: [www.HermanMiller.com/where-to-buy/contact-a-dealer/](http://www.HermanMiller.com/where-to-buy/contact-a-dealer/). Select your market by using the filter.

## 21. Preventive Maintenance

Not applicable

## 22a. Special Attributes.

2022

MillerKnoll earns its 15th consecutive perfect score in Human Rights Campaign Foundation's 2022 Corporate Equality Index.

2021

The OE1 Trolley and Micro Packs are both winners of the 2021 Archiproducts international design competition in the Office Category.

MillerKnoll is named to Michigan's Best and Brightest in Wellness List for 2021.

MillerKnoll is awarded Platinum Rating in CSR by EcoVadis.

OE1 is named winner of Wallpaper's "Most Futuristic Furniture" Category as part of their first-ever Smart Space Awards.

OE1 receives the Workplace category award from Fast Company's 2021 Innovation by Design Awards.

MillerKnoll is named to Fast Company's Annual List of the World's Most Innovative Companies for 2021.

MillerKnoll is recognized as a "Disclosure Leader" by the Chemical Footprint Project (CFP).

# Contract Information

*continued*

MillerKnoll earns its 14th consecutive perfect score in Human Rights Campaign Foundation's 2021 Corporate Equality Index.

MillerKnoll is named to Newsweek's List of America's Most Responsible Companies 2021.

Great Lakes Women's Business Council awards MillerKnoll for "Excellence in Supplier Diversity" in the Best In Class category.

2020

Herman Miller is recognized with a 2020 FSC® Leadership Award for our commitment to responsible forestry management.

Herman Miller is named one of Investor's Business Daily Top 50 Best ESG (Environmental, Social, and Corporate Governance) Companies.

Working Mother names Herman Miller as one of the Best Companies for Dads.

Herman Miller is awarded Platinum Rating in CSR by EcoVadis.

Herman Miller is named Diversity Inc 2020 Noteworthy Company for recruitment efforts, employee development opportunities, leadership accountability, and supplier diversity.

Cosm, by Studio 7.5, receives an iF Gold Award in the Office and Industry category at the iF World Design Awards 2020.

Herman Miller earns its 13th consecutive perfect score in Human Rights Campaign Foundation's 2020 Corporate Equality Index.

2019

Mora System casework, designed by Collective Ten for Herman Miller, receives silver in the Industrial and Life Science Design/ Medical Furniture category at the European Product Design Awards.

Mora System casework is awarded GOOD DESIGN Award 2019, selected for design excellence and innovation.

National Minority Supplier Diversity Council (NMSDC) names Herman Miller as the 2019 Class 1 Corporation of the Year for leadership in supplier diversity.

TIME Magazine names Cosm to its list of 100 Best Inventions of 2019.

Herman Miller is named #12 Overall and #2 in Manufacturing on Investor's Business Daily Best ESG (Environmental, Social, Governance) Companies.

Herman Miller is recognized as a silver level Certified Veteran-Friendly Employer by the Michigan Veterans Affairs Agency for a commitment to veteran hiring, retention, and development.

For the second year in a row, Herman Miller receives a Gold Medal from EcoVadis in recognition of Corporate Social Responsibility achievement measured in categories of environment, labor and human rights, ethics, and sustainable procurement.

Herman Miller is named "Frontrunner" in Chemical Footprint Project by Clean Production Action for our chemical management policies and practices.

Herman Miller receives International Interior Design Association and Contract's best showroom design award in the "Large Showroom" category at NeoCon.

Herman Miller earns its 12th consecutive perfect score in Human Rights Campaign Foundation's 2019 Corporate Equality Index.

Cosm, designed by Studio 7.5 for Herman Miller, receives a Red Dot Best of the Best Product Design Award in the "Office Chairs" category.

Lino, designed by Sam Hecht and Kim Colin for Herman Miller, receives an iF Design Award in the Product Design category.

2018

Overlay, designed by Birsell+Seck, receives Interior Design Best of Year Award in the Partitions and Wall Systems category.

Herman Miller receives the "Corporation of the Year" in consumer products award from the Michigan Minority Supplier Development Council.

Herman Miller achieves the Advanced Excellence in Supplier Diversity Award from the Great Lakes Women's Business Council for doing business with women-owned suppliers.

Great Lakes Women's Business Council awards Herman Miller for "Excellence in Supplier Diversity" in the Advanced category.

WEConnect awards Herman Miller an Honorable Mention for global supplier diversity initiatives.

As one of the top-scoring companies in the industry for the 14th consecutive year, Herman Miller has qualified for inclusion in RobecoSAM's 2018 Sustainability Yearbook and receives the Silver Class distinction for excellent sustainability performance—economic, environmental, and social.

For the 11th consecutive year, Herman Miller receives a perfect score on the Human Rights Campaign Foundation's Corporate Equality Index and is designated as Best Place to Work for LGBT equality. This report evaluates U.S. companies in terms of diversity, particularly LGBT-related policies and practices.

# Contract Information

*continued*

Herman Miller receives Gold Medal from EcoVadis in recognition of Corporate Social Sustainability achievement measured in categories of environment, labor and human rights, ethics, and sustainable procurement. EcoVadis operates the first web-based collaborative platform that allows companies to assess the environmental and social performance of their global suppliers.

Aeron receives Good Design Award from The Chicago Athenaeum Museum of Architecture and Design.

The Cosm Chair wins Mix Interiors "Product of the Year—Task" award.

2017

Herman Miller receives "Rising Star Award" from the Michigan Veterans Affairs Agency (MVAA) for exhibiting new and novel approaches to veteran recruitment and hiring.

Herman Miller is named the 2017 "Corporation of the Year" in the consumer products sector by the Michigan Minority Supplier Development Council (MMSDC).

Herman Miller receives a 2017 SEAL (Sustainability, Environmental Achievement and Leadership) Business Sustainability Award. SEAL Award winners are determined by a holistic methodology measuring applicants against established benchmarks that demonstrate impact and progress toward creating a healthy planet and a sustainable future.

Herman Miller earns the 4 Star Diversity Visionary Award. This award celebrates the diversity, equality, and inclusion work of Herman Miller. The award is presented at the annual Diversity Equity Inclusion Summit, held by the Grand Rapids Area Chamber of Commerce.

Herman Miller is named one of the Best and Brightest in Wellness for the 5th year in a row.

Herman Miller earns its 10th consecutive perfect score in Human Rights Campaign Foundation's 2017 Corporate Equality Index.

As one of the top-scoring companies in the industry for the 13th consecutive year, Herman Miller has qualified for inclusion in RobecoSAM's 2017 Sustainability Yearbook and receives the Bronze Class distinction for excellent economic, environmental, and social sustainability performance. The selection criteria for inclusion in the Yearbook becomes more rigorous each year, and only the top 15 percent of companies within each industry are selected.

Herman Miller once again earns the WorldatWork Seal of Distinction, a unique standard of excellence in work-life effectiveness that shows we provide a distinct, mutually beneficial workplace experience. Herman Miller is the only furniture manufacturer among the 160 companies across the U.S. and Canada honored with the award.

Herman Miller's Mora wins Nightingale Silver Award in the "Furniture Collections" category at the 2017 Healthcare Design Conference.

2016

Herman Miller earns 9th consecutive perfect score in Human Rights Campaign Foundation's Corporate Equality Index and is designated as Best Place to Work for LGBT equality.

Herman Miller is recognized as the Business of the Year at the Annual Meeting of the West Michigan Environmental Action Council (WMEAC) for being an outstanding partner in support of the growth of WMEAC's programs. WMEAC's award goes to the business with the greatest depth and breadth of support, unwavering dedication to environmental action in West Michigan, and a strong corporate commitment to environmental advocacy.

The 101 Best and Brightest Sustainable Companies once again (5th consecutive year) recognizes Herman Miller as an organization in Michigan (and for the first time, in the nation) that exhibits leadership and innovation in its approach to sustainability. This award celebrates those companies that are making their businesses more sustainable, the lives of their employees better, and the community, both locally and globally, more responsible as a whole.

Herman Miller is recognized by the Michigan Minority Supplier Development Council (MMSDC) as a "Corporation of the Year". Additionally, we are the only West Michigan company this year to receive an "Ambassadors Championing Excellence" Award from the MMSDC.

Herman Miller is recognized with "2016 Excellence in Supplier Diversity Award" from the Great Lakes Women's Business Council.

Herman Miller is recognized as a 2016 Healthiest 100 Workplace in America, ranking 49th in the country for our commitment to health and exceptional corporate wellness programming.

Herman Miller receives the Best and Brightest in Wellness Award for the 4th year in a row. This program is an innovative initiative that recognizes and celebrates quality and excellence in health awareness.

The Michigan Occupational Safety and Health Administration renews Herman Miller's "Star" status, the highest workplace safety and health designation, for the Hickory facility in Spring Lake.

Plex Lounge Furniture from Herman Miller receives #MetropolisLikes award at NeoCon 2016.

Plex Lounge Furniture from Herman Miller receives Interior Design HiP Award for Workplace Seating/Lounge at NeoCon 2016.

Herman Miller is named a 2016 Women on Boards Winning Company. Companies selected for this recognition are considered champions of board diversity as 20 percent or more of their board seats are occupied by women.

Herman Miller is selected for inclusion in the 2016 RobecoSAM Sustainability Yearbook for the 12th consecutive year.

Herman Miller again earns the WorldatWork Work-Life 2016 Seal of Distinction. This award identifies organizations that focus on programs that promote work-life balance and overall well-being.

# Contract Information

*continued*

## **22b. Section 508 Compliance**

Not applicable

## **23. Unique Entity ID:** Q2K3MSZ843D8

## **24. Registered in System for Award Management ([www.SAM.gov](http://www.SAM.gov)) Database**

## **25. Cancellation**

Prior to production, no cancellation charge will apply. After production, only actual cost incurred that the contractor can demonstrate if items are sold after 3 months.

## **26. Restocking Charges**

Agencies must notify the contractor for authorization prior to returning any items. The customer agency will be required to pay all packaging and return freight charges. A restocking fee of 50% for Systems Furniture and 35% for all other product lines will be charged for any returned items.

Returns and restocking policy is not applicable to Options products, products under SIN 337127, 339113H, NaughtOne and Nemschoff products.

## **MillerKnoll Inc. Miscellaneous Information:**

Federal ID Number: 38-0837640

Commercial and Government Entity Code (CAGE): 40636

Data Universal Numbering System (DUNS) 00-601-2801

Commercial Entity Code (CEC): 00702290J



<b>Introduction</b>	<b>page 2</b>
<b>Co/Struc® System</b>	<b>3</b>
Structure	5
Work Surfaces	68
Storage	109
Transport	190
<b>Indices</b>	<b>197</b>
By Name	197
By Number	199
<b>Appendices</b>	
20-Day or Less/Assigned Lead-Time Order Information - Healthcare	
Packaging/Keyed-Alike Information-Co/Struc System	
<b>For the latest materials information, refer to the Materials pages on HermanMiller.com.</b>	

# Introduction

## A Note on the Organization of This Book

There are a few signals that will help you understand the organization of this book. Once you know them, you should be able to find your way around easily.

### Like a newspaper, this book is formatted with columns of text and illustrations that run vertically.

Information wraps from one column to the next and continues for as many pages as it needs to.

### Black bars are clues.

A black bar at the top of a page signals the start of information about a product.

This information is divided into two sections, each signaled by a black bar. Product Information includes a description, additional notes, and dimensional drawings. Specification Information, which always begins at the top of a column, includes a series of steps. Each step represents a decision you need to make in order to complete your product number. Steps must be followed sequentially. The number of pages needed to complete information on a product varies; some will be complete in one page; others may fill three or four pages. Just continue going through the steps until you run into the next product, signaled by the black bar at the top of the page.

Black tabs running along the outside edge of the page help you locate what section you are in or what section you want to go to. You can thumb forwards or backwards to find a section; tabs and text appear on both sides of the page.

The additional pages in the back of the book include the appendices and indices. They are signaled by a long black bar running along the outside edge of the page.

### The index is presented two ways.

The first index is in alphabetical order by product name. It can be helpful if you are new to the products and numbering scheme or if you are looking for a type of product, such as work chair or round table.

The second index is by product number. It lists, in alphanumeric order, the base product number, typically a five-digit number.

## General Information

This book is effective October 7, 2024, subject to change without notice. Products may be purchased from authorized Herman Miller dealers who will quote prices upon request. For more information about our products and services or to see a list of dealers, visit [www.HermanMiller.com](http://www.HermanMiller.com).


All prices are list prices. Additional services, such as planning services, design, storage, and installation, are not included and must be added to these prices when the additional services are requested.


Illustrations, specifications, and prices are based on the latest information at the time of publication. The right is reserved to make changes at any time without notice in prices, fabrics, finishes, materials, specifications, and models, and also to discontinue models and finishes.

Dimensions shown in the drawings are engineering drawing dimensions rounded up to the nearest 1/8". Product is shipped set up or knocked down. Products that include "Shipped knocked down" in the description require some assembly. If this statement does not appear in the description, the product is shipped assembled.

Service problems are normally handled by a Herman Miller authorized dealer. Where this cannot be accomplished on a local level, service problems should be referred to the Customer Care Representative for your region, Herman Miller, Inc., Zeeland, Michigan 49464.

### 20-Day or Less/Assigned Lead-Time Shipments

All products and options in this price book not designated by an Assigned Program icon  will ship in **20 business days or less** after being acknowledged by Herman Miller.

Products and options designated by an  are on the assigned lead-time program and may ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

Structure


Work Surfaces

Storage


Transport







Product Information
Description
This assembly attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 25 stand-offs, and 5 spacers.
Notes
To attach assembly to wall, order rail hardware (CO28441, or CO28541) separately.
To finish end, order end cap for standard or seismic rail (CO338) separately.
When seismic components are required, order seismic wall rail (CO479) instead of standard rail (CO282).
Dimensions




Specification Information		
Step 1.		
CO282FF		\$450
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Structure

Rail Hardware

CO284  
CO285




Co/Struc® Structure

Product Information
Description
This package contains enough fasteners to install 12' of standard rail.
Package contains 300 fastener sets.
The 3" sheet metal screws with plastic anchors are #14, flat-head, Phillips-drive screws with 1⅜"-long plastic anchors. The 4" toggle bolts with wall bushings are ¼"-20, flat-head, Phillips-drive toggle bolts.

Specification Information
Step 1.
CO28
Step 2. Hardware
4413" sheet metal screws with plastic anchors
5414" toggle bolts with wall bushings
Prices for Steps 1-2.
CO28441\$336
CO28541\$1044

Adapter Rail


C1610.



**Product Information**

**Description**

This rail attaches to a panel, a module, or wall strips to support Co/Struc® hanging components or a storage case in a straight-line configuration. It does not hang directly on a wall. 2 rails are required to hang a locker from an open panel frame or glazed panel.



**Specification Information**

**Step 1.**

**C1610.**

Step 2. Width	
24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide
60	60" wide


Prices for Steps 1-2.	
C1610. 24	\$300
30	\$312
36	\$328
42	\$356
48	\$378
60	\$422

Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Structure

Equipment Rail, Wall Mounted


C1170.



Co/Struc® Structure

Product Information
Description
This aluminum rail attaches to an architectural wall. It holds work tools or accessories and can be field cut to the appropriate width. The rail includes plastic end caps to accommodate 2 field-cut rail sections.
Notes
Customer must supply attachment hardware for wall conditions.
Order rail-mounted accessories separately:
<ul style="list-style-type: none"><li>• Chart holder (CT730.)</li><li>• Glove box holder (CT726.)</li><li>• Sharps container holder (CT734.)</li><li>• Supply basket (CT720. or CT722.)</li><li>• Suture basket (CT723.)</li><li>• Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)</li><li>• Utility hook (CT776.)</li><li>• Wastebasket and holder (CT728.)</li></ul>
Order work tools separately:
<ul style="list-style-type: none"><li>• Diagonal tray (Y7218.)</li><li>• Mini-shelf (Y7310.)</li><li>• Organizer tray (Y7210.)</li><li>• Paper tray (Y7216.)</li><li>• Rail divider (Y7320.)</li><li>• Vertical tray (Y7217.)</li></ul>
To mount customer-supplied accessories, order rail clamp (CT729.) separately.
Dimensions


Specification Information
Step 1.
C1170.
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
42 42" wide
48 48" wide
60 60" wide
Prices for Steps 1-2.
C1170. 24 \$146
30 \$177
36 \$203
42 \$215
48 \$238
60 \$267
Step 3. End Cap Finish
LU soft white +\$0
WN warm grey neutral +\$0






Equipment Rail, Systems Mounted

C1172.



Product Information
Description
This aluminum rail attaches to a panel, freestanding module, or wall strips and holds work tools or accessories. It includes steel end caps. Attachment hardware is included.
Notes
Order rail-mounted accessories separately:
<ul style="list-style-type: none"><li>• Chart holder (CT730.)</li><li>• Glove box holder (CT726.)</li><li>• Sharps container holder (CT734.)</li><li>• Supply basket (CT720. or CT722.)</li><li>• Suture basket (CT723.)</li><li>• Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)</li><li>• Utility hook (CT776.)</li><li>• Wastebasket and holder (CT728.)</li></ul>
Order work tools separately:
<ul style="list-style-type: none"><li>• Diagonal tray (Y7218.)</li><li>• Mini-shelf (Y7310.)</li><li>• Organizer tray (Y7210.)</li><li>• Paper tray (Y7216.)</li><li>• Rail divider (Y7320.)</li><li>• Vertical tray (Y7217.)</li></ul>
To mount customer-supplied accessories, order rail clamp (CT729.) separately.
Dimensions




Specification Information
Step 1.
C1172.
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
Prices for Steps 1-2.
C1172. 24 \$212
30 \$225
36 \$243
48 \$305
Step 4. End Cap Finish
LU soft white +\$0
WN warm grey neutral +\$0

Co/Struc® Structure

End Cap for Standard or Seismic Rail

CO338




Co/Struc® Structure

Product Information

**Description**  
This end cap provides a finished end for a standard or seismic wall rail and a positive stop for a locker.

**Notes**  
For areas requiring seismic components, order 1 end cap for each end of rail, and order seismic wall rail (CO479) separately.  
For areas not requiring seismic components, order 1 end cap for each end of rail, and order standard rail (CO282) separately.


Dimensions



Specification Information

Step 1.	
CO338FF	\$85

Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Product Information

**Description**  
This aluminum rail attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 26 stand-offs, and 5 spacers.

**Notes**  
Order rail for areas requiring seismic components, when attaching to studs 16" on center, or when attaching to plaster in accordance with OPA-0117 in California.  
To meet the conditions of appropriate OPA, refer to the OPA document for specific hardware. For other applications, order rail hardware (CO28441, or CO28541) separately.  
Order additional seismic components separately:  
• End cap for standard or seismic rail (CO338)  
Field drilling is required for proper installation.


Dimensions



Specification Information

Step 1.		
CO479FF		\$658
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Structure



Co/Struc® Structure

Product Information

Description

This draw rod connects Action Office® Series 1 or 2 panels or Co/Struc® support panels in a straight line. It has a threaded hole in the draw block for attachment to a seismic floor anchor bracket.

Notes

For Action Office Series 1 panels or Co/Struc support panels, specify height of draw rod to match height of panel.

For Action Office Series 2 panels, specify draw rod as follows:

- 32"-high panel, CO215.27
- 39"-high panel, CO215.34
- 47"-high panel, CO215.42
- 53"-high panel, CO215.48
- 62"-high panel, CO215.57
- 67"-high panel, CO215.62
- 85"-high panel, CO215.80

Order floor anchor bracket, Action Office Series 1 (CO481.) or floor anchor bracket, Action Office Series 2 (CO485.) separately.

Dimensions

Specification Information

Step 1.


CO215.

Step 2. Height

27	27" high
34	34" high
42	42" high
48	48" high
57	57" high
62	62" high
80	80" high


Prices for Steps 1-2.

CO215. 27	\$43
34	\$44
42	\$44
48	\$45
57	\$45
62	\$46
80	\$49




Floor Anchor Bracket, Co/Struc® Lab Module

CO568



Product Information
<p>Description</p> <p>This bracket fastens both sides of a module frame to the floor. Each anchor location requires 1 bracket. Finish is light tone. Package contains 10. Hardware to attach module is included. Floor anchor and bolt must be customer supplied.</p>
<p>Notes</p> <p>Order bracket for areas requiring seismic components.</p> <p>Lab modules need to be mounted flush for anchor to install.</p>
<p>Dimensions</p>




Specification Information
<p>Step 1.</p> <p><b>CO568FF</b></p>
<p>\$1455</p>

Co/Struc® Structure

Floor Anchor Bracket, Action Office® Series 1


CO481



Co/Struc® Structure


Product Information
Description
This bracket fastens both sides of an Action Office® Series 1 panel to the floor and is designed for use in areas requiring seismic components. Package contains 10.
Notes
Order the following products separately:
• Draw rod (CO215.)
• Floor anchor adapter kit (CO552.)
• Seismic finished end (A1250.) when using only Action Office Series 2 components
• Seismic hingeable end (AO258.) when using Action Office Series 1 or Co/Struc® components
Customer must supply required bolts and concrete anchors.
Dimensions

Specification Information		
Step 1.		
CO481		\$539
Step 2. Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Floor Anchor Bracket, Action Office® Series 2

CO485



Product Information
Description
This bracket fastens both sides of an Action Office® Series 2 panel to the floor and is designed for use in areas requiring seismic components. Package contains 10.
Notes
Order the following products separately:
• Draw rod (CO215.)
• Floor anchor adapter kit (CO552.)
• Seismic finished end (A1250.)
Customer must supply required concrete anchors.
Dimensions




Specification Information
Step 1.
CO485
\$1193
Step 2. Finish
HF inner tone light
LT light tone
LU soft white
WN warm grey neutral
+\$0
+\$0
+\$0
+\$0

Co/Struc® Structure


Seismic Finished End

A1250.




Co/Struc® Structure

Product Information
Description
This seismic cover finishes the exposed end of a panel. It has a draw block that accepts a seismic bracket to meet the state of California seismic codes, and includes a cable management end cover.
Notes
Order floor anchor bracket, Action Office Series 2 (CO485.) separately.
Dimensions



Specification Information
Step 1.
A1250.
Step 2. Height
39S 39" high
47S 47" high
53S 53" high
62S 62" high
67S 67" high
85S 85" high
Prices for Steps 1-2.
A1250. 39S \$183
47S \$195
53S \$199
62S \$206
67S \$206
85S \$215
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
MT medium tone +\$0
WL warm stone (CP) +\$0
WN warm grey neutral +\$0
Step 4. Cable Management Finish
HF inner tone light +\$0
LU soft white +\$0
MT medium tone +\$0
WL warm stone (CP) +\$0
WN warm grey neutral +\$0






Product Information
<b>Description</b> This kit replaces the bottom block of an Action Office® Series 1 or 2 spacer or 2-way, 3-way, or 4-way connector. It is designed to be used in areas requiring seismic components and has a threaded hole in the draw block for attachment to a seismic floor anchor bracket. Package contains 10.
<b>Notes</b> Order floor anchor bracket, Action Office Series 1 (CO481.) or floor anchor bracket, Action Office Series 2 (CO485.) separately.
<b>Dimensions</b>


Specification Information
<b>Step 1.</b>
<b>CO552</b> <span style="float: right;">\$301</span>

Co/Struc® Structure




Panel Stiffener, Action Office® Series  
1 or Series 2 Top Cap

CO553



Co/Struc® Structure


Product Information
Description
This top cap fits over a standard top cap and spans multiple Action Office® Series 1 or 2 panels. It provides extra rigidity along the panel run. The 12'-long top cap can be field cut to the appropriate length. Attachment hardware is included.
Dimensions



Specification Information		
Step 1.		
CO553FF		\$818
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
MT	medium tone	+\$0
WN	warm grey neutral	+\$0

Wall Strip Cover


CO487.




Product Information
Description
This plastic cover snaps onto a wall strip to cover the slots. It helps with infection control in healthcare environments by preventing bacteria and germs from growing in these slots. It is 84" and must be field cut to match the wall strip height. Package contains 6 covers.
Notes
Order wall strips (AO213. or E1130.) separately.
Dimensions

Specification Information
Step 1.
CO487.
\$655
Step 2. Finish
HF inner tone light
LU soft white
MT medium tone
+\$0
+\$0
+\$0

Co/Struc® Structure





Co/Struc® Structure

Product Information


Description

This open panel has hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

- Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
- To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.
- To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.
- To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.

Dimensions



**Specification Information**

Step 1.

**C1150.**

Step 2. Height

**48** 48" high**80** 80" high

Step 3. Width


**12** 12" wide**24** 24" wide**30** 30" wide**36** 36" wide**48** 48" wide

Prices for Steps 1-3.

	<b>12</b>	<b>24</b>	<b>30</b>	<b>36</b>	<b>48</b>
<b>C1150. 48</b>	\$782	811	872	937	1015
<b>80</b>	\$1076	1145	1205	1278	1393

Step 4. Surface Finish

**HF** inner tone light +\$0**LT** light tone +\$0**LU** soft white +\$0**WN** warm grey neutral +\$0



Co/Struc® Structure

Product Information


Description

This panel has a 48"-high recessed filler and hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

- Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
- To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.
- To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.
- To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.

Dimensions



**Specification Information**

Step 1.

**C1151.**

Step 2. Height

**48** 48" high**80** 80" high

Step 3. Width


**12** 12" wide**24** 24" wide**30** 30" wide**36** 36" wide**48** 48" wide

Prices for Steps 1-3.

	<b>12</b>	<b>24</b>	<b>30</b>	<b>36</b>	<b>48</b>
<b>C1151. 48</b>	\$1038	1288	1426	1549	1741
<b>80</b>	\$1391	1568	1729	1942	2146

Step 4. Surface Finish

**HF** inner tone light +\$0**LT** light tone +\$0**LU** soft white +\$0**WN** warm grey neutral +\$0

**Product Information****Description**

This panel has a vertical raceway for routing and concealing electrical wires and communication cables between the floor and a wire chase rail. It has a 48"-high recessed filler, hanger frames at 1" intervals to support components, and 6 knockouts: 3 at the floor level and 3 at 47" above floor level. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

**Notes**

To supply power from both sides of panel, order wire chase rail (C1313. or C1314.) separately and mount 48" above floor.


For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.


To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.

**Dimensions**



Specification Information				
Step 1.				
<b>C1155.</b>				
Step 2. Height				
<b>48</b>	48" high			
<b>80</b>	80" high			
Step 3. Width				
<b>24</b>	24" wide			
<b>30</b>	30" wide			
<b>48</b>	48" wide			
Prices for Steps 1-3.				
		<b>24</b>	<b>30</b>	<b>48</b>
<b>C1155. 48</b>		\$1826	2023	2407
<b>80</b>		\$2153	2284	2684
Step 4. Surface Finish				
<b>HF</b>	inner tone light			+\$0
<b>LT</b>	light tone			+\$0
<b>LU</b>	soft white			+\$0
<b>WN</b>	warm grey neutral			+\$0

**Product Information****Description**

This panel has a vertical raceway for routing and concealing electrical wires and communication cables between the ceiling and a wire chase rail. It has a 48"-high recessed filler, hanger frames at 1" intervals to support components, and 6 knockouts: 3 at the floor level and 3 at 47" above floor level. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

**Notes**

To carry different types of electrical power (e.g. regular and emergency power) or to carry electrical power and low-voltage cables, vertical raceway can be divided into separate internal raceways; contact local contractor.

To supply power from both sides of panel, order wire chase rail (C1313. or C1314.) separately and mount 48" above floor.


For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.


To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.

To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.

**Dimensions**

Specification Information				
Step 1.				
<b>C1156.</b>				
Step 2. Height				
<b>48</b>	48" high			
<b>80</b>	80" high			
Step 3. Width				
<b>24</b>	24" wide			
<b>30</b>	30" wide			
<b>48</b>	48" wide			
Prices for Steps 1-3.				
		<b>24</b>	<b>30</b>	<b>48</b>
<b>C1156. 48</b>		\$4422	4790	4814
<b>80</b>		\$4842	5056	5237
Step 4. Surface Finish				
<b>HF</b>	inner tone light			+\$0
<b>LT</b>	light tone			+\$0
<b>LU</b>	soft white			+\$0
<b>WN</b>	warm grey neutral			+\$0



Co/Struc® Structure

Product Information




Description

The center point of this 60"-wide panel attaches at a 90° angle to other panels and provides support to a panel run. Each panel has adjustable glides, finished ends, and a 48"-high recessed filler and hanger frames at 1" intervals to support components. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. The panel is available in 2 heights. The 80"-high panel supports a locker; 48"- and 80"-high panels support work surfaces and other 60"-wide components.

Notes


- To support locker, order 2 adapter rails (C1610.) separately.
- To connect to panel at 90° angle, order draw rod (AO215.) separately.
- Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions



Terminal Panel *continued*


Specification Information		
Step 1.		
<b>C1161.</b>		
Step 2. Height		
<b>4860</b>	48" high	
<b>8060</b>	80" high	
Prices for Steps 1-2.		
<b>C1161. 4860</b>		\$2687
<b>8060</b>		\$3303
Step 3. Surface Finish		
<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0



Co/Struc® Structure

Product Information
Description
This connector attaches a support panel to a wall or column for a peninsula configuration.
Notes
Order draw rod (AO215.) separately.
Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.
Dimensions


Specification Information
Step 1.
C1211.
Step 2. Height
4848" high
8080" high
Prices for Steps 1-2.
C1211. 48\$497
80\$594



Step 3. Surface Finish
HFinner tone light+\$0
LTlight tone+\$0
LUsoft white+\$0
WNwarm grey neutral+\$0

2-Way 90° Corner Connector

C1221.



Product Information

Description

This connector joins 2 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.

Notes

Order support panels (C1150., C1151., C1155., C1156.) separately.

To create panel-to-panel connections of same height, order 2 draw rods (AO215.) or 2 seismic draw rods (CO215.) separately.

To create panel-to-panel connection of unequal height or panel-to-corner connection of unequal height, order change-of-height draw rod (AO219.) to match lower height and finished-end change-of-height panel/panel (AO272.) separately.

When connecting panels of unequal heights, specify height of connector to match height of higher panel.

Dimensions

Specification Information

Step 1.

C1221.

Step 2. Height

48	48" high
80	80" high


Prices for Steps 1-2.

C1221. 48	\$497
80	\$595


Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Structure



TOP




SIDE





FRONT

3-Way 90° Corner ConnectorC1231.




Co/Struc® Structure

Product Information
Description
This connector joins 3 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.
Notes
Order support panels (C1150., C1151., C1155., C1156.) separately.
Order 3 draw rods (AO215.) or 3 seismic draw rods (CO215.) separately.
When connecting panels of unequal heights, specify height of connector to match height of higher panel.
Dimensions



Specification Information
Step 1.
C1231.
Step 2. Height
4848" high
8080" high
Prices for Steps 1-2.
C1231. 48\$521
80\$641
Step 3. Surface Finish
HFinner tone light+\$0
LTlight tone+\$0
LUsoft white+\$0
WNwarm grey neutral+\$0







Product Information

**Description**  
This 11"-wide connector joins a module to a support panel. It has a removable cover strip that conceals predrilled holes used for attaching connectors.

**Notes**  
To complete shelf surface at end of module run, order shelf end cap (C1510.E or C1511.E) separately.  
Order panel connector (C1281.) and 3 draw rods (AO215. or CO215.) separately.  
To attach support panel to interface connector (C1261.) of unequal height, order the following products separately:  
• Panel connector (C1281.)  
• Change of height draw rod (AO219.)  
• Change of height panel/panel finished end (AO272.)

Dimensions




Specification Information

**Step 1.**  
**C1261.**

**Step 2. Height**  
**4812** 48" high  
**8012** 80" high

Prices for Steps 1-2.	
<b>C1261. 4812</b>	\$1265
<b>8012</b>	\$1540

Step 3. Surface Finish	
<b>HF</b>	inner tone light +\$0
<b>LT</b>	light tone +\$0
<b>LU</b>	soft white +\$0
<b>WN</b>	warm grey neutral +\$0



Co/Struc® Structure

Product Information

Description


This connector attaches a support panel to an interface connector.

Notes

Specify connector height to match panel height.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions



FRONT

Specification Information

Step 1.

C1281.

Step 2. Height


48	48" high
80	80" high

Prices for Steps 1-2.

C1281.	48	\$117
	80	\$147


Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Product Information
Description
This brace provides additional support to a straight panel run that has components hanging on 1 side. It includes a leveling glide.
Notes
To support 24" deep work surface, specify 21" wide mid-run support; to support 30" deep work surface, specify 27" wide mid-run support. Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
Dimensions

Specification Information		
Step 1.		
C1272.		
Step 2. Height		
48	48" high	
80	80" high	
Step 3. Width		
21	21" wide	
27	27" wide	
Prices for Steps 1-3.		
	21	27
C1272. 48	\$900	996
80	\$1126	1171
Step 4. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0



Mid-Run Support, Double Sided

C1273.



Co/Struc® Structure

Product Information


Description

This brace provides additional support to a straight panel run that has components hanging on both sides. It includes leveling glides.


Notes

To support 24"-deep work surface, specify 41"-wide mid-run support; to support 30"-deep work surface, specify 53"-wide mid-run support. To support 24"-deep work surface on 1 side and 30"-deep work surface on other side, specify 53"-wide mid-run support.  
Order draw rods (AO215.) or seismic draw rods (CO215.) separately.


Dimensions




TOP



FRONT



41" WIDE  
SIDE



53" WIDE  
SIDE

Specification Information

Step 1.

C1273.

Step 2. Height

48	48" high
80	80" high

Step 3. Width

41	41" wide
53	53" wide

Prices for Steps 1-3.


	41	53
C1273. 48	\$985	1145
80	\$1339	1494

Step 4. Surface Finish


HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Support Panel Top Cap


C1274.



Product Information
Description
This 20"-wide aluminum top cap attaches to 2 support panels in a straight connection. It provides additional support for a panel run.
Dimensions




SIDE



FRONT

Specification Information
Step 1.
C1274.20
\$204
Step 2. Surface Finish
HF inner tone light
LT light tone
LU soft white
WN warm grey neutral
+\$0
+\$0
+\$0
+\$0


Co/Struc® Structure




Co/Struc® Structure

Product Information
Description
This 48"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services. The module includes a top shelf and 4 adjustable glides. It has hanger frames at 1" intervals to support components.
Notes
Order access panels separately.
For L, T, or X configuration, order module corner connector separately:
• 2-way connector (C1220.)
• 3-way connector (C1230.)
Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
To attach modules of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.
To enclose the end of a module when a terminal panel (C1122.) is not used, order the end closure assembly (C1451.) separately.
Dimensions

Specification Information
Step 1.
C1112.48
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
Prices for Steps 1-2.
C1112.48 24 \$1893
30 \$2225
36 \$2271
48 \$2476
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0





### Product Information

#### Description

This 80"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services and provides pass-through access. The module includes an interior shelf that is field adjustable to 30", 36", or 48" high. It has hanger frames at 1" intervals to support components and 4 adjustable glides.

#### Notes

Order access panels separately.

For L, T, or X configuration, order module corner connector separately:

- 2-way connector (C1220.)
- 3-way connector (C1230.)


Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach modules of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.

To enclose the top section of a pass-through module, order the end closure panel (C1450.) and 2 change of height panel/panel finished ends (AO272.) separately.

To enclose the end of a module when a terminal panel (C1122., C1124.) is not used, order the end closure assembly (C1451.) separately.

### Dimensions



Co/Struc® Structure

Specification Information

Step 1.

C1114.80

Step 2. Width

- 24 24" wide
- 30 30" wide
- 36 36" wide
- 48 48" wide


Prices for Steps 1-2.

C1114.80	24	\$2687
	30	\$2969
	36	\$3153
	48	\$3346


Step 3. Surface Finish

JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0






Product Information
Description
This 48"-high, hard-surfaced end panel attaches to a freestanding module. It supports a module run or divides space. The panel includes a top shelf end cap, adjustable glides, and finished ends, and comes in 3 widths. The 60"- and 72"-wide panel supports hanging components.
Notes
Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
Dimensions



Specification Information
Step 1.
C1122.48 <span>A</span>
Step 2. Width/Position
36L 36" wide and left-hand panel
36R 36" wide and right-hand panel
60 60" wide
72 72" wide
Prices for Steps 1-2.
C1122.48 36L \$2571
36R \$2571
60 \$3149
72 \$4644
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0

Vertical Chase Cover


C1460.




Co/Struc® Structure


Product Information
Description
This 72"-high tower encloses services between the ceiling and the top shelf of a low module, the interior shelf of a pass-through module, or a 2- or 3-way module corner connector. It includes a ceiling trim plate. The cover cannot be used with a 4-way module corner connector.
Notes
Order module or connector separately:
• 2-way module corner connector (C1220.)
• 3-way module corner connector (C1230.)
• Low module (C1112.)
• Pass-through module (C1114.)
Dimensions

Specification Information		
Step 1.		
C1460.		\$1200
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0






Product Information
Description
This connector starts a peninsula configuration by attaching a low or pass-through module to a wall or column. It includes a shelf end cap.
Notes
For low module, specify 48"-high wall start; for pass-through module, specify 80"-high wall start.
Order draw rods (AO215.) or seismic draw rods (CO215.) separately.
Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.
Dimensions



Specification Information
Step 1.
C1212.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1212. 48 \$847
80 \$956
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0


2-Way Module Corner Connector

C1220.



Co/Struc® Structure


Product Information
Description
This connector joins 2 modules of equal or unequal heights at a 90° angle. It includes a top shelf and filler panels to enclose the utility chase.
Notes
Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach module of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.
Dimensions



Specification Information
Step 1.
C1220.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1220. 48 \$1966
80 \$2359
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0

# 3-Way Module Corner Connector

C1230.



## Product Information

### Description


This connector joins 3 modules of equal or unequal heights at 90° angles. It includes a top shelf and filler panels to enclose the utility chase.

### Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach module of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.

### Dimensions



## Specification Information

### Step 1.

C1230.

### Step 2. Height

48	48" high
80	80" high

### Prices for Steps 1-2.


C1230. 48	\$1966
80	\$2359

### Step 3. Surface Finish

JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0

Module Support Leg


C1275.



Co/Struc® Structure


Product Information
Description
This leg provides additional support for normal loading on a low or pass-through module. It has a leveling glide.
Dimensions

Specification Information		
Step 1.		
C1275.		\$915
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0



Anchor Support Kit


C1276.



Product Information
Description
This kit fastens a module to the floor. It includes a bolt, spacer, and plate.
Package contains 10.
Notes
Anchor support kit is not acceptable for seismic anchoring.
Source anchor component separately.
Dimensions


Specification Information
Step 1.
C1276.
\$452

Co/Struc® Structure



Top Shelf

C1510.




Co/Struc® Structure

Product Information

Description

This 11"-deep shelf encloses the top of a low module service chase and provides a storage shelf. It has a Chemsurf® surface. The shelf can be drilled for fixtures or piping. It cannot attach to a pass-through module.

Dimensions



Specification Information

Step 1.

C1510.12

Step 2. Width

12	12" wide
24	24" wide
30	30" wide
36	36" wide
48	48" wide

Prices for Steps 1-2.

C1510.12	12	\$354
	24	\$475
	30	\$595
	36	\$624
	48	\$640


Step 3. Surface Finish

JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0




Interior Shelf

C1511.



Product Information
Description
This 11"-deep shelf fits inside a pass-through module. It encloses the top of a 48"-high service chase and provides a storage shelf. The shelf has a Chemsurf® surface. It can be drilled for fixtures or piping. Attachment hardware is included.
Notes
To install shelf higher than 48", order interior shelf and support kit (C1515.) separately.
Dimensions




Specification Information
Step 1.
C1511.12
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
Prices for Steps 1-2.
C1511.12 24 \$564
30 \$628
36 \$668
48 \$705
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0

Co/Struc® Structure


Shelf End Cap

C1510.  
C1511.




Co/Struc® Structure

Product Information
Description
This cap finishes the end of an interior or top shelf at the connection of a low or pass-through module and an interface connector. It has a Chemsurf® surface. Attachment hardware is included.
Dimensions




SIDE




TOP


Specification Information
Step 1.
C151
Step 2. Depth/Usage
1.E 7" deep for interior shelf
0.E 12" deep for top shelf
Prices for Steps 1-2.
C1511.E \$314
C1510.E \$307
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0



Product Information
Description
This 11"-deep interior shelf is installed within a pass-through module. Attachment hardware is included.
Dimensions



Specification Information
Step 1.
C1515.12
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
Prices for Steps 1-2.
C1515.12 24 \$1107
30 \$1196
36 \$1330
48 \$1517
Step 3. Surface Finish
JD warm grey neutral Chemsurf® +\$0
KH light tone Chemsurf® +\$0
KR inner tone light Chemsurf® +\$0
LE soft white Chemsurf® +\$0



Product Information

Description

This assembly of 2 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down.

The assembly has 2 30"- or 2 42"-high access panels.


Notes

Order low module (C1112.) or pass-through module (C1114.) separately.

For interior shelf mounted at 36" on pass-through module, specify 30"-high assembly.

For interior shelf mounted at 48" on pass-through module, specify 42"-high assembly.

Dimensions



**Specification Information**

## Step 1.

**C1412.**

## Step 2. Height

**30** 30" high**42** 42" high

## Step 3. Width

**19** 19¼" wide**25** 25¼" wide**31** 31¼" wide**43** 43¼" wide

## Prices for Steps 1-3.


	<b>19</b>	<b>25</b>	<b>31</b>	<b>43</b>
<b>C1412. 30</b>	\$889	1036	1084	1117
<b>42</b>	\$998	1164	1272	1399

## Step 4. Surface Finish

**HF** inner tone light +\$0**LT** light tone +\$0**LU** soft white +\$0**WN** warm grey neutral +\$0

Combination Access Panel Assembly

C1422.



Co/Struc® Structure

Product Information

Description

This assembly of 4 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts (4 per module side) for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down.

The 24"-high assembly has 4 12"-high access panels; the 42"-high assembly has 2 12"- and 2 30"-high access panels.


Notes

Order low module (C1112.) or pass-through module (C1114.) separately.

For interior shelf mounted 30" high on pass-through module, specify 24"-high assembly.

For interior shelf mounted 48" high on pass-through module, specify 42"-high assembly.

Dimensions



# Combination Access Panel Assembly

continued

## Specification Information

Step 1.

**C1422.**

Step 2. Height

**24** 24" high

**42** 42" high

Step 3. Width

**19** 19¼" wide

**25** 25¼" wide

**31** 31¼" wide

**43** 43¼" wide

Prices for Steps 1-3.

	<b>19</b>	<b>25</b>	<b>31</b>	<b>43</b>
<b>C1422. 24</b>	\$1754	1839	1925	1997
<b>42</b>	\$1865	2070	2298	2425

Step 4. Surface Finish

**HF** inner tone light +\$0


**LT** light tone +\$0

**LU** soft white +\$0

**WN** warm grey neutral +\$0

# Low Cup Sink Access Panel Assembly

C1432.



## Product Information

### Description

This assembly attaches to a pass-through module and creates a low cup sink assembly in the lower service chase. It encloses both sides of the service chase and has removable panels for service access. The assembly includes 4 12"-high access panels, 1 shelf support frame, 1 interior shelf, 1 bottom chase cover, and attachment hardware. Each 12"-high panel has 2 knockouts for standard duplex electrical outlets. The assembly has a 15 $\frac{1}{4}$ " space for a cup sink; this divides the service chase into 2 12" segments. Shipped knocked down.


### Notes

Order pass-through module (C1114.) separately.

To enclose service chase of adjacent module, order center closure panel (C1455.1607) separately.

Customer must order sink and fixtures separately and field cut interior shelf to hold sink.

## Dimensions






# Low Cup Sink Access Panel

## Assembly *continued*

Specification Information			
Step 1.			
C1432.42			
Step 2. Width			
19	19¼" wide		
25	25¼" wide		
43	43¼" wide		
Prices for Steps 1-2.			
C1432.42	19		\$2696
	25		\$2944
	43		\$3254
Step 3. Surface Finish			
JD	warm grey neutral Chemsurf®		+\$0
KH	light tone Chemsurf®		+\$0
KR	inner tone light Chemsurf®		+\$0
LE	soft white Chemsurf®		+\$0



End Closure Assembly

C1451.



Co/Struc® Structure

Product Information
Description
This assembly encloses the end of a module when a terminal panel is not used.
The end closure for a low module has a 48"-high panel and 2 48"-high end caps. The full-height end closure for a pass-through module has an 80"-high panel and 2 80"-high end caps.
Notes
Order low module (C1112.) or pass-through module (C1114.) separately.
Dimensions





Specification Information
Step 1.
C145
Step 2. Type
1.4811 end closure for low module
1.8011 full-height end closure for pass-through module
Prices for Steps 1-2.
C1451.4811 \$417
C1451.8011 \$546
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0

Center Closure Panel

C1455.



Product Information
Description
This panel attaches to the inside of a module run and encloses the end of a service chase. Attachment hardware is included.
Notes
To enclose service chase on either end of cup sink opening or 14"-high openings created by partial chases at panel connections, order 16"-high panel.
To enclose lower interior module chase, order 43"-high panel.
To enclose upper section of pass-through module, order 32"-high panel.
Dimensions




Specification Information
Step 1.
C1455.
Step 2. Height
1607 16" high
3207 32" high
4307 43" high
Prices for Steps 1-2.
C1455. 1607 \$939
3207 \$952
4307 \$956
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0

Co/Struc® Structure

## Wire Chase Rail

C1313.  
C1314.



Co/Struc® Structure

### Product Information

#### Description

This aluminum rail attaches to 1 side of an equal-width module or support panel and provides 1 or 2 cable management channels. It cannot be used with an Action Office® Series 1 or 2 connector.

#### Notes


Order connection components separately:

- Wire chase rail bridge (C1326. or C1327.) to connect across mid-run support panel or back side of 3-way connector
- 90° inside corner wire chase rail (C1321. or C1322.) to connect around 90° inside corner


Order wire chase rail covers (C1315., C1316., C1317., or C1318.) and rail end cap (C1331. or C1332.) separately.

Receptacles must be customer supplied.

#### Dimensions



Double Channel



Single Channel

### Specification Information

#### Step 1.

**C131** ☐ A

#### Step 2. Type

3. double channel ☐ A
4. single channel ☐ A

#### Step 3. Width

- 24** 24" wide
- 30** 30" wide
- 36** 36" wide
- 42** 42" wide ☐ A
- 48** 48" wide
- 60** 60" wide
- 72** 72" wide
- 96** 96" wide

#### Prices for Steps 1-3.

<b>C1313.</b>	<b>24</b>	\$174
	<b>30</b>	\$222
	<b>36</b>	\$281
	<b>42</b>	\$328
	<b>48</b>	\$343
	<b>60</b>	\$431
	<b>72</b>	\$512
	<b>96</b>	\$680

<b>C1314.</b>	<b>24</b>	\$92
	<b>30</b>	\$117
	<b>36</b>	\$146
	<b>42</b>	\$164
	<b>48</b>	\$174
	<b>60</b>	\$220
	<b>72</b>	\$258
	<b>96</b>	\$343

#### Step 4. Surface Finish

For double channel (3.) with 24" wide (24)

**AN** clear anodized aluminum +\$0

For double channel (3.) with 30" wide (30)

**AN** clear anodized aluminum +\$0

For double channel (3.) with 36" wide (36)

**AN** clear anodized aluminum +\$0

## Wire Chase Rail *continued*

For double channel (3.) with 42" wide (42)  
AN clear anodized aluminum +\$0

For double channel (3.) with 48" wide (48)  
AN clear anodized aluminum +\$0

For double channel (3.) with 60" wide (60)  
AN clear anodized aluminum +\$0

For double channel (3.) with 72" wide (72)  
AN clear anodized aluminum +\$0

For double channel (3.) with 96" wide (96)  
AN clear anodized aluminum +\$0

For single channel (4.) with 24" wide (24)  
AN clear anodized aluminum +\$0

For single channel (4.) with 30" wide (30)  
AN clear anodized aluminum +\$0

For single channel (4.) with 36" wide (36)  
AN clear anodized aluminum +\$0

For single channel (4.) with 42" wide (42)  
AN clear anodized aluminum +\$0

For single channel (4.) with 48" wide (48)  
AN clear anodized aluminum +\$0


For single channel (4.) with 60" wide (60)  
AN clear anodized aluminum +\$0

For single channel (4.) with 72" wide (72)  
AN clear anodized aluminum +\$0

For single channel (4.) with 96" wide (96)  
AN clear anodized aluminum +\$0

# Wire Chase Rail, 90° Inside Corner

C1321.  
C1322.



Co/Struc® Structure

## Product Information


### Description

This aluminum rail attaches to 2 modules or support panels at a 90° inside corner. It cannot be used with an Action Office® Series 1 or 2 connector.

### Notes

Rail includes 2 6" mitered corner covers. Order wire chase rail covers to cover the remaining length of each rail (C1315., C1316., C1317., or C1318.) separately.

### Dimensions



## Specification Information

### Step 1.

C132 ☐ A

### Step 2. Type

1. double channel ☐ A
2. single channel ☐ A

### Step 3. Width

- 24 24" wide  
30 30" wide  
48 48" wide

### Step 4. Depth

#### For 24" wide (24)

- 24 24" deep  
30 30" deep  
48 48" deep

#### For 30" wide (30)

- 30 30" deep  
48 48" deep

#### For 48" wide (48)

- 48 48" deep

### Prices for Steps 1-4.

	24	30	48
C1321. 24	\$596	641	765
30	—	\$679	804
48	—	—	\$927
C1322. 24	\$401	431	492
30	—	\$451	512
48	—	—	\$576

### Step 5. Surface Finish

#### For double channel (1.) with 30" wide (30) with 48" deep (48)

AN clear anodized aluminum +\$0

#### For double channel (1.) with 48" wide (48) with 48" deep (48)

AN clear anodized aluminum +\$0

#### For single channel (2.) with 24" wide (24) with 24" deep (24)

AN clear anodized aluminum +\$0

## Wire Chase Rail, 90° Inside Corner

*continued*

---

*For single channel (2.) with 24" wide (24) with 30" deep (30)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For single channel (2.) with 24" wide (24) with 48" deep (48)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For single channel (2.) with 30" wide (30) with 30" deep (30)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For single channel (2.) with 30" wide (30) with 48" deep (48)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For single channel (2.) with 48" wide (48) with 48" deep (48)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For double channel (1.) with 24" wide (24) with 24" deep (24)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For double channel (1.) with 24" wide (24) with 30" deep (30)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---

*For double channel (1.) with 24" wide (24) with 48" deep (48)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---


*For double channel (1.) with 30" wide (30) with 30" deep (30)*

AN	clear anodized aluminum	+\$0
----	-------------------------	------

---


Wire Chase Rail Bridge

C1326.  
C1327.



Co/Struc® Structure

Product Information
Description
This bridge joins wire chase rails across a mid-run support panel or on the back side of a 3-way 90° connector. It allows the continuation of a cable management channel and includes a blank cover. The bridge cannot be used with an Action Office® Series 1 or 2 connector.
Notes
Order 1 of the following products separately: <ul style="list-style-type: none"><li>• 3-way 90° connector (C1231.)</li><li>• Mid-run support (C1272. or C1273.)</li></ul>
Dimensions




Specification Information
Step 1.
C132 <input type="checkbox"/> A
Step 2. Type
6. double channel <input type="checkbox"/> A
7. single channel <input type="checkbox"/> A
Prices for Steps 1-2.
C1326. \$154
C1327. \$154
Step 3. Surface Finish
AN clear anodized aluminum +\$0




Wire Chase Rail End Cap

C1331.  
C1332.



Product Information
Description
This end cap finishes the exposed end of a wire chase rail. Package contains 2.
Dimensions




Specification Information
Step 1.
C133 <span>A</span>
Step 2. Type
1. double channel <span>A</span>
2. single channel <span>A</span>
Prices for Steps 1-2.
C1331. \$138
C1332. \$138
Step 3. Surface Finish
AN clear anodized aluminum +\$0

Co/Struc® Structure

Wire Chase Rail Cover

C1315.  
C1316.  
C1317.  
C1318.




Co/Struc® Structure

Product Information

Description  
This cover encloses a wire chase rail.

Notes  
Specify type to cover wire chase rail's power and/or data outlets.  
Specify covers in 6"- or 12"-wide increments to match rail width.

Dimensions




Specification Information

Step 1.  
C131

Step 2. Type  
5.06 6"-wide plain cover   
5.12 12"-wide plain cover   
6. duplex cutout cover   
7. GFI cutout cover   
8. gang cover

Prices for Steps 1-2.	
C1315.06	\$35
C1315.12	\$44
C1316.	\$47
C1317.	\$49
C1318.	\$49


Step 3. Surface Finish  
AN clear anodized aluminum +\$0




Product Information
Description
This 1½"-wide cover wraps around the outside of a wire chase channel to cover gaps between adjoining channels and covers.
Notes
Order wire chase and covers separately.
Dimensions

Specification Information
Step 1.
C1358. <span>A</span>
Step 2. Type
D double channel
T triple channel
Prices for Steps 1-2.
C1358. D \$44
T \$62

Co/Struc® Structure



**Product Information****Description**

This work surface hangs from a panel, module, or wall strips and is available as laminate or solid surface. The work surface, when supported by work surface supports has a top-loading capacity of up to 200 pounds. The solid surface option, when supported by sink products has a top-loading capacity of 250 pounds. Appropriate sink products and work surface supports must be specified and ordered separately.

**Notes**

For laminate option (L), order standard duty work surface supports (C2140.S) separately.

For solid surface option (K), order 1 of the following products separately:

- ADA sink support, systems mounted (C2520.)
- Sink module (C2512.)
- Sink unit for storage (C2513.)

- Work surface support, standard duty (C2140.S)

Work surface with solid surface option (K) will not accept a storage case (C5040.) without storage case leg (C5050.).


To support a 30" deep work surface installed on wall strips, a storage case (C5040.) with storage case leg (C5050.) is required. Order separately.

Work surface with solid surface option (K) has a matte gloss level and includes a backsplash.

When hanging work surface with laminate option (L) from wall strips, order recommended backsplash (C2910.) separately.

Order optional products separately:

- Drawer bearer package (CO228)
- Gripper package (CO266), for laminate (L) option only
- Keyboard/Mouse Support, Adjustable (Y7710.)
- Lockable drawer bearer (CO396)

**Dimensions****Specification Information****Step 1.****C2142.****Step 2. Depth**

<b>24</b>	24" deep
<b>30</b>	30" deep

**Step 3. Width**

<b>24</b>	24" wide
<b>30</b>	30" wide
<b>36</b>	36" wide
<b>48</b>	48" wide
<b>60</b>	60" wide
<b>72</b>	72" wide
<b>96</b>	96" wide

**Step 4. Material**

For 24" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), or 60" wide (60)

<b>L</b>	high - pressure laminate
<b>K</b>	solid surface

For 72" wide (72) or 96" wide (96)

<b>L</b>	high - pressure laminate
----------	--------------------------

**Prices for Steps 1-4.**

			<b>L</b>	<b>K</b>
<b>C2142.</b>	<b>24</b>	<b>24</b>	\$167	2046
		<b>30</b>	\$206	2375
		<b>36</b>	\$277	2979
		<b>48</b>	\$369	3370
		<b>60</b>	\$535	4034
		<b>72</b>	\$671	—
		<b>96</b>	\$859	—
	<b>30</b>	<b>24</b>	\$177	2291
		<b>30</b>	\$256	2659
		<b>36</b>	\$321	3331
		<b>48</b>	\$512	3755
		<b>60</b>	\$725	4489
		<b>72</b>	\$891	—
		<b>96</b>	\$1065	—

## Standard-Duty Work Surface *continued*


### Step 5. Finish

#### *For high - pressure laminate (L)*

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

#### *For solid surface (K)*

<b>CQ</b>	glacier white	+\$0
<b>QB</b>	bisque	+\$0
<b>1S</b>	savannah	+\$52
<b>1U</b>	sahara	+\$52
<b>4I</b>	sandstone	+\$52
<b>58</b>	cameo white	+\$50
<b>ANC</b>	antarctica	+\$50
<b>OQ</b>	silt	+\$52
<b>SVG</b>	silver gray	+\$50
<b>CRT</b>	concrete	+\$200
<b>DVC</b>	dove	+\$200



### Product Information

#### Description

This work surface is supported by heavy duty work surface supports, allowing it to hang from a panel, module, or wall strips. It has a top-loading capacity of up to 400 pounds. It can hold drawers and includes drawer bearers. Appropriate heavy duty work surface supports must be specified and ordered separately.

The work surface can be specified as laminate, Chemsurf®, stainless steel, Corian®, or phenolic resin. A 51"-wide surface has a 3" overhang on the left or right side to fill the gap created by a mid-run support, C-series 3-way connector, or interface connector. A 96"-wide surface includes 1 additional support bracket. The 24"-wide surface cannot accept components.

#### Notes

Order heavy duty work surface supports (C2140.H, C2140.W, C214.R) separately.

For stainless steel option (S), order 1 of the following products separately:

- Heavy duty work surface support, with component adapter (C2140.W)
- Heavy duty work surface support, component adapter ready (C2140.R)

Work surface with phenolic resin option (P) or Corian option (K) must use work surface supports (C2140.) heavy duty, with component adaptor (W) or heavy duty, component adaptor ready (R) options.

Finish color on phenolic resin option (P) is black.

Work surface with solid surface option (K) has a matte gloss level.

Work surface with phenolic resin option (P) or Corian option (K) includes a backsplash.

Work surface with phenolic resin option (P) is 1" thick and includes spacers for transitioning to adjacent, non-resin surfaces.

To support a 30" deep work surface installed on wall strips, a storage case (C5040.) with storage case leg (C5050.) is required. Order separately.

When hanging a work surface with laminate option (L) or Chemsurf® option (C) from wall strips, order recommended backsplash (C2910.) separately.

Work surface supports the following quantities of drawers, C frames, or storage cabinets:

- 24" wide: 1 drawer
- 30", 36" wide: 1 drawer, 1 C frame or 1 storage cabinet
- 48", 51", and 60" wide: 2 drawers, 2 C frames, or 2 storage cabinets
- 72" wide: 3 drawers, 3 C frames, or 3 storage cabinets
- 96" wide: 4 drawers, 4 C frames, or 4 storage cabinets


24"-wide work surface will not accept a storage case (C5040.) or storage unit (C5010.) or C frame (CO205).

Work surface with phenolic resin option (P) or stainless steel option (S) will not accept a keyboard/mouse support, adjustable (Y7710.).

Order optional products separately:

- Drawer (CO207, CO208, or CO209)
- Keyboard tray (Y7710.), except for stainless steel (S) and phenolic resin (P) surfaces
- Pencil drawer (A0480.)
- Storage case (C5040.) and adapter rail (C1610.), except for phenolic resin (P) surfaces
- Storage unit (C5010.) and component adaptor (C2193.)

### Dimensions



# Heavy-Duty Work Surface *continued*

## Specification Information

### Step 1.

**C2141.**

### Step 2. Depth

**24** 24" deep

**30** 30" deep

### Step 3. Width

**24** 24" wide

**30** 30" wide

**36** 36" wide

**48** 48" wide

**5L** 51" wide, left overhang

**5R** 51" wide, right overhang

**60** 60" wide

**72** 72" wide

**96** 96" wide

### Step 4. Material

*For 24" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide (60), or 72" wide (72)*

**L** high-pressure laminate

**C** Chemsurf®

**S** stainless steel ☐

**P** phenolic resin

**K** corian

*For 51" wide, left overhang (5L) or 51" wide, right overhang (5R)*

**L** high-pressure laminate

**C** Chemsurf®

**S** stainless steel ☐

**P** phenolic resin

**K** corian

*For 96" wide (96)*

**L** high-pressure laminate

**C** Chemsurf®

**S** stainless steel ☐

## Prices for Steps 1-4.

		L	C	S	P	K
<b>C2141.</b>	<b>24 24</b>	\$227	305	3057	1084	1797
	<b>30</b>	\$242	357	3289	1283	2008
	<b>36</b>	\$314	453	3926	1476	2447
	<b>48</b>	\$428	572	4691	1867	3299
	<b>5L</b>	\$511	660	5109	2066	3746
	<b>5R</b>	\$511	660	5109	2066	3746
	<b>60</b>	\$594	735	6113	2162	3748
	<b>72</b>	\$735	895	6975	2460	4629
	<b>96</b>	\$936	1118	7809	—	—
	<b>30 24</b>	\$293	346	3441	1182	2047
	<b>30</b>	\$305	465	3882	1476	2295
	<b>36</b>	\$415	619	4935	1676	2776
	<b>48</b>	\$557	786	5809	2162	3691
	<b>5L</b>	\$581	812	6281	2460	4200
	<b>5R</b>	\$581	812	6281	2460	4200
	<b>60</b>	\$766	995	7455	2555	4192
	<b>72</b>	\$936	1188	8586	2948	5204
	<b>96</b>	\$1149	1451	9167	—	—

## Step 5. Finish

*For high-pressure laminate (L)*

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

*For Chemsurf® (C)*

<b>JD</b>	warm grey neutral Chemsurf®	+\$0
<b>KH</b>	light tone Chemsurf®	+\$0
<b>KR</b>	inner tone light Chemsurf®	+\$0
<b>LE</b>	soft white Chemsurf®	+\$0

## Heavy-Duty Work Surface *continued*

Co/Struc® Work Surfaces


*For corian (K)*

<b>CQ</b>	glacier white	+\$0
<b>QB</b>	bisque	+\$0
<b>1S</b>	savannah	+\$380
<b>1U</b>	sahara	+\$380
<b>4I</b>	sandstone	+\$380
<b>58</b>	cameo white	+\$370
<b>ANC</b>	antarctica	+\$370
<b>IQ</b>	aurora	+\$380
<b>OQ</b>	silt	+\$380
<b>SVG</b>	silver gray	+\$100
<b>CRT</b>	concrete	+\$740
<b>DVC</b>	dove	+\$740



# Work Surface Support

C2140.



## Product Information

### Description

This set of supports allows a standard- or heavy duty work surface to hang from a panel, module, or wall strips.

### Notes

Supports for standard duty surfaces include the following:

Surface Width—Supports

24"-60" wide—2 supports (1 left and 1 right)

72" wide—3 supports (1 left, 1 middle, and 1 right)

96" wide—4 supports (1 left, 2 middle, and 1 right)

Supports for heavy duty surfaces include the following:

Surface Width—Supports


24"-72" wide—2 supports (1 left and 1 right)

96" wide—3 supports (1 left, 1 middle, and 1 right)

For standard duty option (S), order standard duty work surface (C2142.) separately.

For heavy duty options (W or R), order appropriate heavy duty work surface (C2141.) separately.

## Dimensions



## Specification Information

### Step 1.

**C2140.**

### Step 2. Type

<b>S</b>	standard duty
<b>W</b>	heavy duty, with component adaptor
<b>R</b>	heavy duty, component adaptor ready

### Step 3. Depth

<b>24</b>	for 24"-deep work surface
<b>30</b>	for 30"-deep work surface

### Step 4. Width

<b>24</b>	24" wide
<b>30</b>	30" wide
<b>36</b>	36" wide
<b>48</b>	48" wide
<b>60</b>	60" wide
<b>72</b>	72" wide
<b>96</b>	96" wide


### Prices for Steps 1-4.

		24	30	36	48	60	72
<b>C2140. S</b>	<b>24</b>	\$252	252	252	252	252	378
	<b>30</b>	\$252	252	252	252	252	378
	<b>W 24</b>	\$634	674	708	746	814	886
	<b>30</b>	\$693	739	770	808	876	948
	<b>R 24</b>	\$456	456	456	456	456	502
	<b>30</b>	\$518	518	518	518	518	570
							<b>96</b>
<b>C2140. S</b>	<b>24</b>						\$502
	<b>30</b>						\$502
	<b>W 24</b>						\$1222
	<b>30</b>						\$1317
	<b>R 24</b>						\$747
	<b>30</b>						\$849

### Step 5. Finish

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

Rail-Hanging Work SurfaceC2410.




Product Information

**Description**  
This 24"-deep rectangular work surface hangs on a standard rail. It has predrilled holes for hanging drawers, suspended pedestals, or keyboard trays. The work surface can be installed, relocated, or removed without tools and without disturbing other rail components. The support arm finish matches the surface finish.

**Notes**  
Order optional work surface backsplash (C2910.) separately.  
Order optional keyboard/mouse support, adjustable (Y7710.) separately.

Dimensions




Specification Information

Step 1.		
C2410.24		
Step 2. Width		
48	48" wide	
72	72" wide	
Prices for Steps 1-2.		
C2410.24	48	\$1643
	72	\$2032
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Pass-Through Cover


C2399.




Product Information
Description
This cover provides a continuous surface for 2 back-to-back squared-edge work surfaces hung on an open panel frame. It is 48" wide but can be field cut to fit a narrower opening.
Dimensions

Specification Information		
Step 1.		
C2399.48		\$151
Step 2. Finish		
LU	soft white	+\$0

Co/Struc® Work Surfaces






Co/Struc® Work Surfaces

Product Information
Description
This assembly attaches to a support panel, a module, or wall strips and provides support for a C frame and components hanging on the C frame.
Attachment hardware is included. Shipped knocked down.
The assembly holds the following quantities of C frames:
• 24" and 30" wide: 1 C frame
• 48" and 60" wide: 2 C frames
• 72" wide: 3 C frames
Notes
For 51"-wide work surface, order 48"-wide C frame support.
Dimensions


Specification Information
Step 1.
C2192.
Step 2. Work Surface Depth
24 for 24"-deep work surface
30 for 30"-deep work surface
Step 3. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
60 60" wide
72 72" wide

Prices for Steps 1-3.							
		24	30	36	48	60	72
C2192.	24	\$828	851	872	912	919	935
	30	\$912	971	1015	1033	1043	1092
Step 4. Surface Finish							
HF	inner tone light						+\$0
LT	light tone						+\$0
LU	soft white						+\$0
WN	warm grey neutral						+\$0




Component Adapter Kit

C2193.




Product Information
Description
This kit attaches to a heavy-duty work surface support and allows C frames to hang below the work surface. Attachment hardware is included.
The kit holds the following quantities of C frames or storage units:
• 30" wide: 1 C frame or storage unit
• 48" and 60" wide: 2 C frames or storage units
• 72" wide: 3 C frames or storage units
• 96" wide: 4 C frames or storage units
Notes
The kit attaches to a work surface support (C2140.) with heavy duty, component adapter ready option (R).
For 51"-wide work surface, order 48"-wide kit.
96"-wide kit is shipped in 2 pieces.
Order storage components separately:
• C frame (CO205)
• Storage unit (C5010.)
Dimensions




Specification Information
Step 1.
C2193.
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
60 60" wide
72 72" wide
96 96" wide
Prices for Steps 1-2.
C2193. 24 \$239
30 \$286
48 \$338
60 \$409
72 \$504
96 \$642
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0

Work Surface Backsplash

C2910.




Product Information
Description
This plastic barrier attaches to the back edge of a C-series rectangular or corner work surface or to a process table and creates a backsplash. Package contains 6.
Dimensions



Specification Information
Step 1.
C2910.
Step 2. Width
72 72" wide
Prices for Steps 1-2.
C2910. 72 \$624
Step 3. Surface Finish
For 72" wide (72)
BK black +\$0
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0

Backsplash Filler


C2920.




Product Information
Description
This component attaches to a Casework work surface. It fills the gap created when Casework base cabinets are placed against terminal panels or wall strips. Package contains 6. Attachment hardware is included.
Notes
Order work surface (CG211. or CG212.) or corner work surface (CG218. or CG219.) separately.
Dimensions

Specification Information		
Step 1.		
C2920.96		\$226
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Work Surfaces



Heavy-Duty Corner Work SurfaceC2215.



Product Information

Description

This 48"-wide work surface hangs from modules, C-series support panels, or wall strips and provides additional top-loading capacity. It includes 1 drawer bearer package. The work surface cannot attach to Action Office® Series 1 or Series 2 connectors. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top. The support arm finish matches the surface finish.

Notes


When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard/mouse support, adjustable (Y7710.) separately.


Dimensions



TOP



TOP



FRONT



# Heavy-Duty Corner Work Surface

continued

## Specification Information

### Step 1.

**C2215.** ☐

### Step 2. Depth

**2448** 24" deep ☐

**3048** 30" deep ☐

### Prices for Steps 1-2.

<b>C2215. 2448</b>	\$1126
<b>3048</b>	\$1305

### Step 3. Surface Finish

#### Laminate Top

##### For 24" deep (2448)

<b>HF</b>	inner tone light <input type="checkbox"/>	+\$0
<b>LT</b>	light tone <input type="checkbox"/>	+\$0
<b>LU</b>	soft white <input type="checkbox"/>	+\$0
<b>WN</b>	warm grey neutral <input type="checkbox"/>	+\$0

#### Chemsurf® Top

##### For 24" deep (2448)

<b>JD</b>	warm grey neutral Chemsurf® <input type="checkbox"/>	+\$174
<b>KH</b>	light tone Chemsurf® <input type="checkbox"/>	+\$174
<b>KR</b>	inner tone light Chemsurf® <input type="checkbox"/>	+\$174
<b>LE</b>	soft white Chemsurf® <input type="checkbox"/>	+\$174

#### Laminate Top

##### For 30" deep (3048)


<b>HF</b>	inner tone light <input type="checkbox"/>	+\$0
<b>LT</b>	light tone <input type="checkbox"/>	+\$0
<b>LU</b>	soft white <input type="checkbox"/>	+\$0
<b>WN</b>	warm grey neutral <input type="checkbox"/>	+\$0

#### Chemsurf® Top

##### For 30" deep (3048)

<b>JD</b>	warm grey neutral Chemsurf® <input type="checkbox"/>	+\$199
<b>KH</b>	light tone Chemsurf® <input type="checkbox"/>	+\$199
<b>KR</b>	inner tone light Chemsurf® <input type="checkbox"/>	+\$199
<b>LE</b>	soft white Chemsurf® <input type="checkbox"/>	+\$199

# Heavy-Duty Corner Work Surface, C2216. Action Office® Series 1 or 2 90° Corner Connector



## Product Information

### Description

This 48"-wide work surface hangs from panels or wall strips connected by an Action Office® Series 1 or 2 90° corner connector. It provides additional top-loading capacity and includes 1 drawer bearer package. The support arm finish matches the surface finish. The work surface cannot attach to a C-series 90° corner connector. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top.

### Notes

When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard/mouse support, adjustable (Y7710.) separately.

### Dimensions



## Specification Information

### Step 1.

**C2216.** ☐

### Step 2. Depth

**2448** 24" deep ☐

**3048** 30" deep ☐

### Prices for Steps 1-2.

<b>C2216. 2448</b>	\$1145
<b>3048</b>	\$1326

### Step 3. Surface Finish

#### Laminate Top

##### For 24" deep (2448)

<b>HF</b>	inner tone light <input type="checkbox"/>	+\$0
<b>LT</b>	light tone <input type="checkbox"/>	+\$0
<b>LU</b>	soft white <input type="checkbox"/>	+\$0
<b>WN</b>	warm grey neutral <input type="checkbox"/>	+\$0

#### Chemsurf® Top

##### For 24" deep (2448)

<b>JD</b>	warm grey neutral Chemsurf® <input type="checkbox"/>	+\$174
<b>KH</b>	light tone Chemsurf® <input type="checkbox"/>	+\$174
<b>KR</b>	inner tone light Chemsurf® <input type="checkbox"/>	+\$174
<b>LE</b>	soft white Chemsurf® <input type="checkbox"/>	+\$174

#### Laminate Top

##### For 30" deep (3048)

<b>HF</b>	inner tone light <input type="checkbox"/>	+\$0
<b>LT</b>	light tone <input type="checkbox"/>	+\$0
<b>LU</b>	soft white <input type="checkbox"/>	+\$0
<b>WN</b>	warm grey neutral <input type="checkbox"/>	+\$0


#### Chemsurf® Top

##### For 30" deep (3048)

<b>JD</b>	warm grey neutral Chemsurf® <input type="checkbox"/>	+\$199
<b>KH</b>	light tone Chemsurf® <input type="checkbox"/>	+\$199
<b>KR</b>	inner tone light Chemsurf® <input type="checkbox"/>	+\$199
<b>LE</b>	soft white Chemsurf® <input type="checkbox"/>	+\$199

# Process Table, Electric Height Adjustment

C4313.



## Product Information

### Description

This height-adjustable process table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from 26½" to 41½". The table base is equipped with a motor, height-control box, and 8½' power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 3" locking casters. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

The table includes 1 stretcher to support a C frame or storage unit.

### Notes


Order storage components separately:

- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Reagent shelf (C4326.)
- Storage unit (C5010.)

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.

### Dimensions



## Specification Information

### Step 1.

C4313.

### Step 2. Depth

24 24" deep

30 30" deep

### Step 3. Width

48 48" wide

60 60" wide

72 72" wide

### Step 4. Surface Material

L high-pressure laminate

C Chemsurf® A

### Prices for Steps 1-4.

		L	C
C4313.	24 48	\$8213	8661
	60	\$8482	8928
	72	\$8839	9287
	30 48	\$8368	8812
	60	\$8542	8988
	72	\$8894	9340

### Step 5. Casters/Glides

GA glides +\$0

CA casters +\$318

### Step 6. Top/Frame Finish

#### For high-pressure laminate (L)

HF inner tone light +\$0

LT light tone +\$0

LU soft white +\$0

WN warm grey neutral +\$0

#### For Chemsurf® (C)

JD warm grey neutral Chemsurf® +\$0


KH light tone Chemsurf® +\$0

KR inner tone light Chemsurf® +\$0


LE soft white Chemsurf® +\$0

Reagent Shelf


C4326.  
C4327.



Product Information
Description
This metal shelf attaches to a workstation or process table to provide storage for reagents below the surface.
Notes
Order workstation (C4310. or 4311.) separately.
Order process table (C4313.) separately.
Dimensions



For Electric Process Table




For Workstation

Specification Information			
Step 1.			
C432			
Step 2. Type			
6.	for electric process table		
7.	for workstation		
Step 3. Depth			
For electric process table (6.)			
16	16" deep		
For workstation (7.)			
18	18" deep		
Step 4. Width			
48	48" wide		
60	60" wide		
72	72" wide		
Prices for Steps 1-4.			
	48	60	72
C4326. 16	\$1206	1254	1338
C4327. 18	\$1174	1304	1315
Step 5. Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0
LU	soft white		+\$0
WN	warm grey neutral		+\$0

## Workstation, Electric Height Adjustment

C4310.



### Product Information

#### Description

This 30"-deep height-adjustable table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from 26½" to 41½". The table base is equipped with a motor, height-control box, and 8½' power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

A riser frame attaches to the table base to provide upright support for hanging components. The table includes 1 stretcher to support a C frame or storage unit. The 48"-wide workstation has 1 set of drawer bearers; 60"- and 72"-wide workstations have 2 sets of drawer bearers.

#### Notes

60"- and 72"-wide riser frames require a riser frame divider for overhead storage components. See Co/Struc Planning Guide for information.

Order storage components (to hang on stretcher) separately:

- B-style shelf (C0521.)
- C frame (C0205)
- Drawer (C0207, C0208, or C0209)

Order workstation storage components separately:

- Reagent shelf (C4327.)

Order overhead components (to hang on riser frame) separately:


- B-style shelf (C3011.)
- Equipment rail (C1172.)
- Flat panel mount (C4324.)
- Marker board (C4323.)
- Power strip (C4322.)
- Riser frame divider (C4320.)
- Riser frame end cover (C4321.)
- Shelf unit (C5020. or C5021.)

Order Action Office® overhead components separately. See Co/Struc Planning Guide for information.

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.

### Dimensions



# Workstation, Electric Height Adjustment *continued*

Co/Struc® Work Surfaces

## Specification Information

### Step 1.

**C4310.30**

### Step 2. Width

**48** 48" wide  
**60** 60" wide  
**72** 72" wide

### Step 3. Surface Material

**L** high-pressure laminate  
**C** Chemsurf®

### Prices for Steps 1-3.

		<b>L</b>	<b>C</b>
<b>C4310.30</b>	<b>48</b>	\$9211	9641
	<b>60</b>	\$9333	9765
	<b>72</b>	\$9758	10190

### Step 4. Casters/Glides

**GA** glides +\$0  
**CA** casters +\$380

### Step 5. Top/Frame Finish

#### *For high-pressure laminate (L)*


**HF** inner tone light +\$0  
**LT** light tone +\$0  
**LU** soft white +\$0  
**WN** warm grey neutral +\$0

#### *For Chemsurf® (C)*

**JD** warm grey neutral Chemsurf® +\$0  
**KH** light tone Chemsurf® +\$0  
**KR** inner tone light Chemsurf® +\$0  
**LE** soft white Chemsurf® +\$0

# Workstation, Manual Height Adjustment

C4311.



## Product Information

### Description

This 30"-deep height-adjustable work station has a laminate or Chemsurf® top. It is pin adjustable in 1" increments and has a 15" travel range from 26½" to 41½". The table has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

A riser frame attaches to the table base to provide upright support for hanging components. The table includes 1 stretcher to support a C frame or storage unit. The table includes 1 set of drawer bearers; 60"- and 72"-wide workstations have 2 sets of drawer bearers.

### Notes

60"- and 72"-wide riser frames require a riser frame divider for overhead storage components. See Co/Struc Planning Guide for information.

Order storage components (to hang on stretcher) separately:

- B-style shelf (C0521.)
- C frame (C0205)
- Drawer (C0207, C0208, or C0209)
- Storage unit (C5010.)

Order workstation storage components separately:

- Reagent shelf (C4327.)

Order overhead components (to hang on riser frame) separately:


- Equipment rail (C1172.)
- Flat panel mount (C4324.)
- Marker board (C4323.)
- Power strip (C4322.)
- Riser frame divider (C4320.)
- Riser frame end cover (C4321.)

Order Action Office® overhead components separately. See Co/Struc Planning Guide for information.

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.

## Dimensions



# Workstation, Manual Height Adjustment *continued*

Co/Struc® Work Surfaces

## Specification Information

### Step 1.

**C4311.30**

### Step 2. Width

**48** 48" wide  
**60** 60" wide  
**72** 72" wide

### Step 3. Surface Material

**L** high-pressure laminate  
**C** Chemsurf®

### Prices for Steps 1-3.

		<b>L</b>	<b>C</b>
<b>C4311.30</b>	<b>48</b>	\$6843	7299
	<b>60</b>	\$7000	7453
	<b>72</b>	\$7324	7916

### Step 4. Casters/Glides

**GA** glides +\$0  
**CA** casters +\$380

### Step 5. Top/Frame Finish

#### For high-pressure laminate (L)

**HF** inner tone light +\$0  
**LT** light tone +\$0  
**LU** soft white +\$0  
**WN** warm grey neutral +\$0


#### For Chemsurf® (C)

**JD** warm grey neutral Chemsurf® +\$0  
**KH** light tone Chemsurf® +\$0  
**KR** inner tone light Chemsurf® +\$0  
**LE** soft white Chemsurf® +\$0



# Work Table, Electric Height Adjustment

C4312.



## Product Information

### Description

This 2-sided, 42"-deep height-adjustable table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from 26½" to 41½". The table base is equipped with a motor, height-control box, and 8½' power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

The table includes 2 stretchers (one on each side) to support C frames and/or storage units. 48"- and 60"-wide tables have 4 sets of drawer bearers.

### Notes




Order storage components separately:

- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (C5010.)

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.

### Dimensions



## Specification Information

### Step 1.

**C4312.42**

### Step 2. Width

<b>48</b>	48" wide
<b>60</b>	60" wide

### Step 3. Surface Material

<b>L</b>	high-pressure laminate
<b>C</b>	Chemsurf®

### Prices for Steps 1-3.

		<b>L</b>	<b>C</b>
<b>C4312.42</b>	<b>48</b>	\$9511	10012
	<b>60</b>	\$9812	10310

### Step 4. Casters/Glides

<b>GA</b>	glides	+\$0
<b>CA</b>	casters	+\$380


### Step 5. Top/Frame Finish

#### For high-pressure laminate (L)

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

#### For Chemsurf® (C)

<b>JD</b>	warm grey neutral Chemsurf®	+\$0
<b>KH</b>	light tone Chemsurf®	+\$0
<b>KR</b>	inner tone light Chemsurf®	+\$0
<b>LE</b>	soft white Chemsurf®	+\$0



### Product Information

#### Description

This freestanding table base is available at seated height (29" high) and standing height (36" high), has 2 C-legs, a steel modesty panel, and 1" adjustable glides. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

#### Notes

If casters are desired at initial purchase, Options is recommended as a solution to maintain 29" and 36" heights. Casters are also available as a service part if replacing glides with casters is desired in the future.

At lowest glide adjustment, actual height of table with glides and work surface is 28 $\frac{3}{4}$ " for a 29" seated-height table base and 35 $\frac{3}{4}$ " for a 36" standing-height table base. Actual height of table with casters and work surface is 31.4" for a 29" seated-height table base and 38.4" for a 36" standing-height table base. Order casters as a service part.

Clearance from front of table to modesty panel is 24 $\frac{1}{4}$ " for a 30"-deep work surface and 18 $\frac{1}{4}$ " for a 24"-deep work surface.

Order work table surface (C4700.) separately:

- 24" and 30" depths, and 48", 51" (left and right), and 60" widths are compatible with table bases 24" deep x 48" wide.
- 60" wide work table surface overhangs on center of a 48"-wide table base.
- 51" wide (left and right) work table surfaces overhang left or right of a 48"-wide table base.
- Order drawer bearers separately, with work table surface.

Order separately:





- Work table reagent shelf (C4900.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).

To accommodate a refrigerator, a 36"-high table base is recommended.

Refrigerators 21.34" wide or less can be used with 1 set of adjacent drawer bearers for a 3" or 6" drawer (CO207, CO208) only.


Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.

### Dimensions



Basic Work Table Base *continued*

Specification Information			
Step 1.			
<b>C460</b>	<input type="text" value="A"/>		
Step 2. Glides			
<b>G.</b>	glides	<input type="text" value="A"/>	
Step 3. Depth			
<b>24</b>	24" deep	<input type="text" value="A"/>	
Step 4. Width			
<b>48</b>	48" wide	<input type="text" value="A"/>	
Step 5. Height			
<b>29</b>	29" high	<input type="text" value="A"/>	
<b>36</b>	36" high	<input type="text" value="A"/>	
Prices for Steps 1-5.			
		<b>29</b>	<b>36</b>
<b>C460G. 24 48</b>		\$1544	1630
Step 6. Base/Modesty Panel Finish			
<b>HF</b>	inner tone light	<input type="text" value="A"/>	+\$0
<b>LT</b>	light tone	<input type="text" value="A"/>	+\$0
<b>LU</b>	soft white	<input type="text" value="A"/>	+\$0
<b>WN</b>	warm grey neutral	<input type="text" value="A"/>	+\$0



### Product Information

#### Description

This freestanding table base is available at seated height (29" high) and standing height (36" high) and has stretchers to support C frames, drawers, and storage units. Table base comes with 1" adjustable glides. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

#### Notes

If casters are desired at initial purchase, Options is recommended as a solution to maintain 29" and 36" heights. Casters are also available as a service part if replacing glides with casters is desired in the future.

At lowest glide adjustment, actual height of table with glides and work surface is 28 $\frac{3}{4}$ " for a 29" seated-height table base and 35 $\frac{3}{4}$ " for a 36" standing-height table base. Actual height of table with casters and work surface is 31.4" for a 29" seated-height table base and 38.4" for a 36" standing-height table base.

Order work table surface (C4700.) separately:

- 24" deep × 36" wide base is only compatible with a 24" deep × 36" wide work table surface.
- 24" deep × 48" wide base is compatible with a 24" deep × 48", 51" (left and right), 60" wide work table surface.
- 24" deep × 72" wide base is only compatible with a 24" deep × 72" wide work table surface.
- 30" deep × 36" wide base is only compatible with a 30" deep × 36" wide work table surface.
- 30" deep × 48" wide base is compatible with a 30" deep × 48", 51" (left and right), 60" wide work table surface.
- 30" deep × 48S" wide base is only compatible with a 36" deep × 72" wide work table surface.
- 30" deep × 72" wide base is only compatible with a 30" deep × 72" wide work table surface.

60" wide work table surface (C4700.) overhangs on center of a 48"-wide table base.

51" wide (left and right) work table surfaces (C4700.) overhang left or right of a 48"-wide table base.

Order drawer bearers separately, with work table surface (C4700.).

Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.

Order storage components separately:

- Drawer (CO207, CO208, CO209).
- Storage unit (C5010.).
- Work table reagent shelf (C4900.).
- C frame (CO205).

Maximum C frames or storage units using maximum number of stretchers and a 36"-high table base:

- For 36"-wide base: 2 C frames or 1 storage unit.
- For 48"-wide base: 4 C frames or 2 storage units.
- For 48S"-wide × 30"-deep base: 8 C frames (4 front/4 back) or 4 storage units.


- For 72"-wide base: 6 C frames or 3 storage units.

Order technology support and cord management accessories separately:

- Connect™ - S300 (Y1323.xxxxS).
- Flo® single-screen monitor arm (Y91171.CM).
- Tablet support (Y91173.).
- Laptop mount (Y91177.).
- LE series keyboard solutions (Y7710.1).
- JAW™ adjustable CPU holder assembly (Y91103.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).
- Cable trough (Y1200.).

Order Y1200.30 for a C461G.3036; order Y1200.42 for a C461G.3048.

## Dimensions



### Specification Information

Step 1.

C461 

A
---

### Step 2. Glides

**G.** glides A

### Step 3. Depth

24 24" deep

30 30" deep

### Step 4. Width

For 24" deep (24)

36 36" wide A

48      48" wide   ☐ A

72 72" wide A

For 30" deep (30)

36 36" wide A

48      48" wide   ☐ A

72 72" wide ☐ A

**48S** 48" wide, symmetrical support A

### Step 5. Height

29 29" high

36 36" high

### Step 6. Number of Stretchers

*For 29" high (29)*

A upper placement only A

For 36" high (36)

A upper placement only A

**B** upper and lower placement A

### Prices for Steps 1-6.


		29A	36A	36B	
C461G.	24	36	\$1772	1879	2183
		48	\$1791	1897	2199
		72	\$1940	2047	2485
30	36		\$1815	1916	2222
		48	\$1836	1935	2242
		72	\$1985	2089	2527
		48S	\$2223	2406	2783

Co/Struc® Work Surfaces

Step 7. Base/Stretcher Finish			
HF	inner tone light	A	+\$0
LT	light tone	A	+\$0
LU	soft white	A	+\$0
WN	warm grey neutral	A	+\$0

## Adjustable Work Table Base

C462C.  
C462G.



### Product Information

#### Description

This freestanding 1" pin-adjustable table base has stretchers to support C frames, drawers, and storage units. Table base comes with 1" adjustable glides or 4" locking casters. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

#### Notes

Table with glides and work surface adjusts from 27" high to 36" high. Table with casters adjusts from 30" high to 38" high.

Order work table surface (C4700.) separately:

- 24" deep × 36" wide base is only compatible with a 24" deep × 36" wide work table surface.
- 24" deep × 48" wide base is compatible with a 24" deep × 48", 51" (left and right), 60" wide work table surface.
- 24" deep × 72" wide base is only compatible with a 24" deep × 72" wide work table surface.
- 30" deep × 36" wide base is only compatible with a 30" deep × 36" wide work table surface.
- 30" deep × 48" wide base is compatible with a 30" deep × 48", 51" (left and right), 60" wide work table surface.
- 30" deep × 48S" wide base is only compatible with a 36" deep × 72" wide work table surface.
- 30" deep × 72" wide base is only compatible with a 30" deep × 72" wide work table surface.

60" wide work table surface (C4700.) overhangs on center of a 48"-wide table base.

51" wide (left and right) work table surfaces (C4700.) overhang left or right of a 48"-wide table base.

Order drawer bearers separately, with work table surface (C4700.).

Order storage components separately:

- Drawer (CO207, CO208, CO209).
- Storage unit (C5010.).
- Work table reagent shelf (C4900.).
- C frame (CO205).

Maximum C frames or storage units using maximum number of stretchers and a 36"-high table base:

- For 36"-wide base: 2 C frames or 1 storage unit.
- For 48"-wide base: 4 C frames or 2 storage units.
- For 48S"-wide × 30"-deep base: 8 C frames (4 front/4 back) or 4 storage units.
- For 72"-wide base: 6 C frames or 3 storage units.

For tables adjusted below 32" high, C frames hang from top stretcher only, and storage units cannot mount.


Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.

Order technology support and cord management accessories separately:

- Connect™ - S300 (Y1323.xxxxS).
- Flo® single-screen monitor arm (Y91171.CM).
- Tablet support (Y91173.).
- Laptop mount (Y91177.).
- LE series keyboard solutions (Y7710.1).
- JAW™ adjustable CPU holder assembly (Y91103.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).
- Cable trough (Y1200.).

Order Y1200.30 for a C462x.3036; order Y1200.42 for a C462x.3048.

### Dimensions



# Adjustable Work Table Base *continued*

## Specification Information

### Step 1.

**C462** ☐

### Step 2. Casters/Glides

**G.** glides ☐

**C.** casters ☐

### Step 3. Depth

**24** 24" deep ☐

**30** 30" deep ☐

### Step 4. Width

#### For 24" deep (24)

**36** 36" wide ☐

**48** 48" wide ☐

**72** 72" wide ☐

#### For 30" deep (30)

**36** 36" wide ☐

**48** 48" wide ☐

**72** 72" wide ☐

**48S** 48" wide, symmetrical support ☐

### Step 5. Number of Stretchers

**A** upper placement only ☐

**B** upper and lower placement ☐

### Prices for Steps 1-5.

	A	B
<b>C462G. 24 36</b>	\$2756	3063
<b>48</b>	\$2770	3084
<b>72</b>	\$2929	3376
<b>30 36</b>	\$2782	3086
<b>48</b>	\$2800	3110
<b>72</b>	\$2954	3401
<b>48S</b>	\$3201	3596

	A	B
<b>C462C. 24 36</b>	\$2864	3175
<b>48</b>	\$2880	3192
<b>72</b>	\$3039	3486
<b>30 36</b>	\$2887	3194
<b>48</b>	\$2907	3217
<b>72</b>	\$3062	3510
<b>48S</b>	\$3307	3704

### Step 6. Base Finish


**HF** inner tone light ☐ +\$0

**LT** light tone ☐ +\$0

**LU** soft white ☐ +\$0

**WN** warm grey neutral ☐ +\$0





### Product Information

#### Description

This work table surface attaches to work table bases (C460G., C461G., C462G., C462C.). Predrilled holes accept drawer bearers ordered with the work surface. Surface material is available in laminate, Chemsurf®, Corian®, or phenolic resin.

#### Notes

Drawer bearers must be field installed.

Only used on work table bases C460G., C461G., C462G., C462C., ordered separately. See table base information for work surface compatibilities and top load weight capacities.


Work table surface supports the following maximum pairs of drawer bearers (each pair has 1 left and 1 right bearer):

- For 36"-wide table: 1 pair.
- For 48"-wide table: 2 pair.
- For 48S"-wide table: 4 pair (2 pair on either side).
- For 60"-wide table: 2 pair.
- For 72"-wide table: 3 pair.

Phenolic resin is 1" thick. When ordering adjustable work table base (C462C.) with phenolic resin top (C4700.xxxxP) that needs to be level with tables (C462C.) using laminate, Chemsurf, or Corian tops (C4700.L, C4700.K, C4700.C), order phenolic shim service part.

See Co/Struc® System planning guide for capacity information.

### Dimensions



# Work Table Surface, Rectangular

continued

## Specification Information

Step 1.

**C4700.** ☐

Step 2. Depth

**24** 24" deep ☐

**30** 30" deep ☐

**36** 36" deep ☐

Step 3. Width

For 24" deep (24)

**36** 36" wide ☐

**48** 48" wide ☐

**51L** 51" wide, left overhang ☐

**51R** 51" wide, right overhang ☐

**60** 60" wide ☐

**72** 72" wide ☐

For 30" deep (30)

**36** 36" wide ☐

**48** 48" wide ☐

**51L** 51" wide, left overhang ☐

**51R** 51" wide, right overhang ☐

**60** 60" wide ☐

**72** 72" wide ☐

For 36" deep (36)

**72** 72" wide ☐

Step 4. Surface Material

**L** high-pressure laminate ☐

**C** Chemsurf® ☐

**K** corian ☐

**P** phenolic resin ☐

Prices for Steps 1-4.

	L	C	K	P
<b>C4700. 24 36</b>	\$350	468	1180	1467
<b>48</b>	\$381	510	1472	1947
<b>51L</b>	\$480	600	1769	2361
<b>51R</b>	\$480	628	1769	2361
<b>60</b>	\$492	612	1770	2363
<b>72</b>	\$535	751	2071	2848

<b>30 36</b>	\$391	502	1389	1806
<b>48</b>	\$441	546	1730	2397
<b>51L</b>	\$540	683	2082	2909
<b>51R</b>	\$540	683	2082	2909
<b>60</b>	\$535	749	2072	2900
<b>72</b>	\$603	840	2419	3494
<b>36 72</b>	\$797	1241	2824	5334

Step 5. Top/Edge Finish

For high-pressure laminate (L)

**HF** inner tone light ☐ +\$0

**LT** light tone ☐ +\$0

**LU** soft white ☐ +\$0

**WN** warm grey neutral ☐ +\$0

For Chemsurf® (C)

**JD** warm grey neutral Chemsurf® ☐ +\$0

**KH** light tone Chemsurf® ☐ +\$0

**KR** inner tone light Chemsurf® ☐ +\$0

**LE** soft white Chemsurf® ☐ +\$0

For corian (K)

**CQ** glacier white ☐ +\$0

**QB** bisque ☐ +\$0

**1S** savannah ☐ +\$410

**1U** sahara ☐ +\$410

**4I** sandstone ☐ +\$410

**58** cameo white ☐ +\$400

**ANC** antarctica ☐ +\$400

**IQ** aurora ☐ +\$410

**OQ** silt ☐ +\$410

**SVG** silver gray ☐ +\$100

**CRT** concrete ☐ +\$525

**DVC** dove ☐ +\$525

For phenolic resin (P)

**BK** black ☐ +\$0

Step 6. Pairs of Drawer Bearers

For 36" wide (36)

**NB** without drawer bearers ☐ +\$0

**1D** 1 pair drawer bearers ☐ +\$16

## Work Table Surface, Rectangular

*continued*

*For 48" wide (48), 51" wide, left overhang (51L), 51" wide, right overhang (51R), or 60" wide (60)*

<b>NB</b>	without drawer bearers <input type="checkbox"/>	+\$0
<b>1D</b>	1 pair drawer bearers <input type="checkbox"/>	+\$16
<b>2D</b>	2 pair drawer bearers <input type="checkbox"/>	+\$31

*For 24" deep (24) or 30" deep (30) with 72" wide (72)*


<b>NB</b>	without drawer bearers <input type="checkbox"/>	+\$0
<b>1D</b>	1 pair drawer bearers <input type="checkbox"/>	+\$16
<b>2D</b>	2 pair drawer bearers <input type="checkbox"/>	+\$31
<b>3D</b>	3 pair drawer bearers <input type="checkbox"/>	+\$47

*For 36" deep (36) with 72" wide (72)*

<b>NB</b>	without drawer bearers <input type="checkbox"/>	+\$0
<b>1D</b>	1 pair drawer bearers <input type="checkbox"/>	+\$16
<b>2D</b>	2 pair drawer bearers <input type="checkbox"/>	+\$31
<b>3D</b>	3 pair drawer bearers <input type="checkbox"/>	+\$47
<b>4D</b>	4 pair drawer bearers <input type="checkbox"/>	+\$62

Work Table Reagent Shelf

C4900.



Co/Struc® Work Surfaces

Product Information

Description

This high-pressure laminate shelf attaches to work table bases to provide storage below the surface.

Notes

Compatible with basic work table base (C460G.), standard work table base (C461G.xxxxxxxA), and adjustable work table base with upper stretcher placement only (C462x.xxxxxxA).


Compatible only with 36"-high table bases, due to installation clearance needed.

16" deep shelf (C4900.16) is compatible with 24"- and 30"-deep table bases.

18" deep shelf (C4900.18) is compatible with only 30"-deep table bases.

Order warm grey neutral (WN) reagent shelf finish to complement a metallic champagne (CN) table base. Order cool grey neutral (CL) reagent shelf finish to complement a metallic silver (MS) table base.

Reagent shelf (C4900.) rests on table base legs and is not height adjustable.



Specification Information

Step 1.

C4900.

Step 2. Depth

16	16" deep	<input type="text"/>
18	18" deep	<input type="text"/>

Step 3. Width

36	36" wide	<input type="text"/>
48	48" wide	<input type="text"/>
72	72" wide	<input type="text"/>

Step 4. Surface Material

L	high-pressure laminate	<input type="text"/>
---	------------------------	----------------------

Prices for Steps 1-4.


			L
C4900.	16	36	\$594
		48	\$669
		72	\$966
	18	36	\$623
		48	\$696
		72	\$999

Step 5. Top/Edge Finish

HF	inner tone light	<input type="text"/>	+\$0
LT	light tone	<input type="text"/>	+\$0
LU	soft white	<input type="text"/>	+\$0
WN	warm grey neutral	<input type="text"/>	+\$0
CL	cool grey neutral	<input type="text"/>	+\$0

Riser Frame Divider


C4320.




Product Information
Description
This piece divides a workstation riser frame into smaller sections and provides upright riser support for hanging components. The dividers bolt to the riser frame at predrilled attachment locations. Finish is black.
Notes
Order workstation (C4310. or C4311.) separately.
Dimensions

Specification Information
Step 1.
C4320. \$510

Co/Struc® Work Surfaces



Riser Frame End CoverC4321.



Product Information

Description

This metal cover finishes the exposed hardware on the end of a riser frame.

Notes

Specify 2 end covers per riser frame.  
Specify end cover finish to match the corresponding finish of workstation (C4310. or C4311.) to which cover will be attached.  
Order workstation (C4310. or C4311.) separately.

Dimensions


Specification Information


Step 1.

C4321.	\$135
--------	-------

Step 2. Finish


HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0






Product Information
Description
This power strip mounts horizontally on hanger clips and provides single-circuit, 15-amp, 100-volt outlets on a riser frame. It has a lighted on/off switch with circuit breaker and an 8' cord. It is UL listed. Mounting hardware is included.
Power strips have the following number of outlets:
Width—Outlets
24"—4 outlets
30"—5 outlets
36"—6 outlets
48"—8 outlets
60"—10 outlets
72"—12 outlets
Notes
To mount power strip inside riser frame (vertically), specify 36" width.
Specify power strip finish to match the corresponding finish of workstation (C4310. or C4311.) to which it will be attached.
Order workstation (C4310. or C4311.) separately.
For lab modules (C1112., C1114.) and support panels (C1150., C1151.), power strip can be specified as part of the product.
Dimensions

Specification Information
Step 1.
C4322.
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
48 48" wide
60 60" wide
72 72" wide
Prices for Steps 1-2.
C4322. 24 \$588
30 \$632
36 \$677
48 \$699
60 \$788
72 \$832
Step 3. Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0



Marker Board

C4323.



Product Information


**Description**

This 2-sided white board mounts between 2 riser frame dividers or a riser frame and a riser frame divider. Both sides are magnetic with white erasable writing surfaces. Includes a black-texture marker tray on one side only. Mounting hardware is included.

**Notes**

Order workstation (C4310. or C4311.) separately.

Dimensions



Specification Information

Step 1.

C4323.

Step 2. Width

24	24" wide
30	30" wide
36	36" wide


Prices for Steps 1-2.

C4323.	24	\$560
	30	\$620
	36	\$660



Flat Panel Mount


C4324.




Product Information
<p><b>Description</b></p> <p>This arm mounts to a workstation riser frame and holds a flat panel display monitor. The arm meets VESA standards (Video Electronic Standards Association - Mounting Interface Standards) with a 100 mm or 75 mm mounting pattern. It provides 360° monitor rotation and 60° horizontal and vertical tilt range. Finish is black. Mounting hardware is included. Finish is black.</p>
<p><b>Notes</b></p> <p>Order workstation (C4310. or C4311.) separately.</p>
<p><b>Dimensions</b></p>

Specification Information
<p>Step 1.</p> <p><b>C4324.</b></p>
<p>\$845</p>

Co/Struc® Work Surfaces





Co/Struc® Work Surfaces



Product Information

Description  
This component hangs from a module, a support panel, or wall strips and holds a 12"-high sink. It includes a front panel and side panels to conceal the bottom of the sink and left and right work surface supports. Attachment hardware is included. Shipped knocked down.

- The module accepts the following cutout dimensions:
- Depth:
- 24": 18 1/4"-deep maximum cutout
  - 30": 24 1/4"-deep maximum cutout
- Width:
- 24": 19 5/8"-wide maximum cutout
  - 30": 25 5/8"-wide maximum cutout
  - 48": 43 5/8"-wide maximum cutout
  - 60": 55 5/8"-wide maximum cutout

Notes  
Sink must be customer supplied and field installed.  
Order standard duty work surface (C2142.) with solid surface option (K) separately.

Dimensions



Specification Information

Step 1.  
C2512.

Step 2. Depth	
24	24" deep
30	30" deep


Step 3. Width	
24	24" wide
30	30" wide
36	36" wide
48	48" wide
60	60" wide

Prices for Steps 1-3.					
	24	30	36	48	60
C2512. 24	\$836	973	1022	1070	1164
30	\$994	1164	1226	1281	1435


Step 4. Surface Finish		
91	white (CP)	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
76	light brown walnut	+\$0

Sink Unit for Storage

C2513.




Product Information
Description
This 24"-deep, 48"-wide structure hangs from 2 24"-wide panels or freestanding modules or from 3 wall strips. It supports drawer storage on the left or right side and a sink unit on the other side. The maximum cutout for the sink is 18 1/4" deep x 19 5/8" wide. The surface includes a front panel and side panels to conceal the bottom of the sink and left and right work surface supports. Attachment hardware is included. Shipped knocked down.
Notes
Order storage components for open space opposite sink unit separately:
• Drawer bearers (CO228 or CO396)
• Gripper package (CO266) and C frame (CO205) for Co/Struc® drawers
• L cart (CO212) and single counter top (CO555) for C frames and drawers when surface is mounted 38" above floor
• Storage case (C5040.2422) and 24"-wide Co/Struc adapter rail (C1610.24)
Sink must be customer supplied and field installed.
Storage unit (C5010.) cannot mount under sink unit.
Order standard duty work surface (C2142.) with solid surface option (K) separately.
Dimensions



Specification Information
Step 1.
C2513.24
\$1079
Step 2. Surface Finish
HF inner tone light
LT light tone
LU soft white
WN warm grey neutral
91 white (CP)
+\$0
+\$0
+\$0
+\$0
+\$77

Co/Struc® Work Surfaces

ADA Sink Support, Systems Mounted C2520.



Product Information

Description

This sink support hangs from a module or wall strips. It supports a 5½"-deep sink and a solid surface or alternative material top. The sink support has a 4"-high laminate front trim and a removable front panel to allow access to plumbing. Attachment hardware is included. Shipped knocked down.

Notes

Order standard duty work surface (C2142.) with solid surface option (K) separately.  
Sink must be customer supplied and field installed. A Just SL-ADA-1921-A-GR sink, with drain located in rear center, fits in a 24" x 24" sink support.

Dimensions

Specification Information

Step 1.

C2520.

Step 2. Depth

24	24" deep
30	30" deep

Step 3. Width



24	24" wide
30	30" wide
48	48" wide
60	60" wide

Prices for Steps 1-3.

	24	30	48	60
C2520. 24	\$1417	1595	1770	1945
30	\$1505	1680	1945	2123


Step 4. Surface Finish

91	white (CP)	+\$123
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Storage Unit

C5010.




Product Information
Description
This cabinet hangs below a heavy-duty work surface or process table. It has 2 side panels, a back panel, and a top and bottom panel.
Notes
Order heavy-duty work surface (C2141.), process table (C4313.), standard work table base (C461G.), adjustable lab work table base (C462G., C462C.), work table (C4312.), or workstation (C4310. or C4311.) separately.
Order shelves (C5015.1622) separately; shelves adjust in 1¼" increments.
Order A-size drawers (CO207) or B-size drawers (CO208) and lockable drawer bearer (C5048.) separately.
Order storage unit door (C5016. or C5017.) separately.
Dimensions



Specification Information		
Step 1.		
C5010.1622		\$1000
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
91	white (CP)	+\$62

Co/Struc® Storage

Storage Unit ShelfC5015.



Product Information

Description

This shelf mounts inside a storage unit and adjusts in 1¼" increments.

Notes

Order storage unit (C5010.) separately.

Dimensions


Specification Information

Step 1.

C5015.1622	\$280
------------	-------

Step 2. Surface Finish


HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
91	white (CP)	+\$16



Co/Struc® Storage

## Storage Unit Door

C5016.  
C5017.



### Product Information

#### Description

This door attaches to the front of a storage unit. Attachment hardware is included.


#### Notes

Co/Struc® drawers cannot be specified behind doors.

When specifying 19"- or 22"-high lockable door that has drawers above it, specify interior filler (C5009.) to prevent access to cabinet below.

For general information on locks, see Keyed-Alike Information in Appendices.

#### Dimensions



### Specification Information

#### Step 1.

**C501**

#### Step 2. Hinge Location

- 6. left-hand hinges
- 7. right-hand hinges

#### Step 3. Height

- 1922** 19" high
- 2222** 22" high
- 2522** 25" high

#### Step 4. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

#### Prices for Steps 1-4.


		<b>N</b>	<b>L</b>	<b>K</b>
<b>C5016.</b>	<b>1922</b>	\$361	451	451
	<b>2222</b>	\$414	503	503
	<b>2522</b>	\$482	569	569
<b>C5017.</b>	<b>1922</b>	\$361	451	451
	<b>2222</b>	\$414	503	503
	<b>2522</b>	\$482	569	569

#### Step 5. Surface Finish

<b>HF</b>	inner tone light			+\$0
<b>LT</b>	light tone			+\$0
<b>LU</b>	soft white			+\$0
<b>WN</b>	warm grey neutral			+\$0
<b>91</b>	white (CP)			+\$31

#### Step 6. Pull Option/Finish


<b>NN</b>	none			+\$0
<b>BK</b>	black			+\$0
<b>NH</b>	brushed nickel			+\$0



Co/Struc® Storage

Product Information
Description
This filler mounts inside a storage unit and provides closure to the cabinet below when Co/Struc drawers are installed. It mounts at either 19" or 22" high.
Notes
Order storage unit (C5010.) separately.
Filler is used with A-size drawers (CO207) or B-size drawers (CO208).
Dimensions


Specification Information		
Step 1.		
C5009.1622		\$280
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
91	white (CP)	+\$16





Lockable Drawer Bearer


C5048.



Product Information
Description
These drawer bearers attach Co/Struc® A-, B-, or C-size drawers in a storage unit. They have integral stops to prevent drawers from falling when pulled out. Mounting hardware is included. Finish is black. Package contains 1 pair.
Notes
Order A-size drawers (CO207), B-size drawers (CO208), or C-size drawers (CO209) separately.
Can be used with locking or non-locking drawers.
Dimensions


Specification Information
Step 1.
C5048.
\$122

Co/Struc® Storage



Storage Case

C5040.



Product Information

**Description**

This case mounts under a heavy-duty rectangular or corner work surface. It has 2 side panels, a back panel, and a top and bottom panel. Interior components can be height adjusted in 1¼" increments.

Notes

Order heavy-duty work surface (C2141.) or heavy-duty corner work surface (C2215. or C2216.) separately.


To support case mounted under rectangular work surface, order adapter rail (C1610.)

Order optional products separately:

- Storage case bottle drawer (C5045.)
- Storage case door (C5052. or C5053.)
- Storage case drawer (C5043. or C5044.)
- Storage case pullout shelf (C5042.)
- Storage case shelf (C5041.)

For top slot of case, specify 16"-deep drawer or shelf to clear support beam under heavy-duty work surface.

Dimensions



Specification Information

Step 1.

**C5040.**

Step 2. Depth

<b>2422</b>	23" deep
<b>3022</b>	29" deep

Prices for Steps 1-2.


<b>C5040. 2422</b>	\$1097
<b>3022</b>	\$1205

Step 3. Surface Finish

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0
<b>91</b>	white (CP)	+\$72

Storage Case Shelf


C5041.



**Product Information**  
Description  
This shelf mounts inside a storage case and adjusts in 1¼" increments.

**Notes**  
To clear support beam on heavy-duty work surface, shelf must be positioned at least 2½" below the surface.

**Dimensions**



**Specification Information**  
Step 1.  
C5041.


**Step 2. Depth**  
2422 23" deep  
3022 29" deep

**Prices for Steps 1-2.**  
C5041. 2422 \$346  
3022 \$364

**Step 3. Surface Finish**  
HF inner tone light +\$0  
LT light tone +\$0  
LU soft white +\$0  
WN warm grey neutral +\$0  
91 white (CP) +\$21


Co/Struc® Storage

Storage Case Pullout ShelfC5042.



Product Information
Description
This pullout shelf is 21¼" deep and mounts inside a storage case. It adjusts in 1¼" increments.
Notes
To clear support beam on heavy-duty work surface, shelf must be positioned at least 2½" below the surface.
Dimensions


Specification Information
Step 1.
C5042.2422\$485
Step 2. Surface Finish
HFinner tone light+\$0
LTlight tone+\$0
LUsoft white+\$0
WNwarm grey neutral+\$0
91white (CP)+\$33



## Storage Case Drawer

C5043.

C5044.



## Product Information


### Description

This drawer mounts inside a storage case and is available with or without a lock. The 22"-deep drawer cannot be used in the top position of a storage case.

## Notes

For general information on locks, see Keyed-Alike Information in Appendices.

## Dimensions



### Specification Information

Step 1.

## C504

### Step 2. Depth

- 16" deep
- 22" deep

### Step 3. Height

For 16" deep (3.)

- 0422 small, 5" high  
0822 medium, 7½" high

For 22" deep (4.)

- 0422 small, 5" high  
0822 medium, 7½" high  
1222 large, 10" high

### Step 4. Lock

- |   |                   |
|---|-------------------|
| N | no lock           |
| L | keyed differently |
| K | keyed alike       |

## Prices for Steps 1-4.

	N	L	K
C5043. 0422	\$459	758	758
0822	\$576	864	864
C5044. 0422	\$508	669	669
0822	\$607	899	899
1222	\$639	933	933

### Step 5. Surface Finish


- |           |                   |       |
|-----------|-------------------|-------|
| <b>HF</b> | inner tone light  | +\$0  |
| <b>LT</b> | light tone        | +\$0  |
| <b>LU</b> | soft white        | +\$0  |
| <b>WN</b> | warm grey neutral | +\$0  |
| <b>91</b> | white (CP)        | +\$41 |

### Step 6. Pull Option/Finish

- |           |                |      |
|-----------|----------------|------|
| <b>NN</b> | none           | +\$0 |
| <b>BK</b> | black          | +\$0 |
| <b>NH</b> | brushed nickel | +\$0 |

Storage Case Bottle Drawer

C5045.




Product Information

Description

This pullout drawer mounts in a storage case and has a partial-height front for easy access to contents. It includes 3 dividers.

Dimensions



Specification Information

Step 1.

C5045.

Step 2. Height

0822 5" high

1222 7 1/2" high

Prices for Steps 1-2.

C5045. 0822 \$1014

1222 \$1080

Step 3. Surface Finish

HF inner tone light +\$0

LT light tone +\$0

LU soft white +\$0

WN warm grey neutral +\$0

76 light brown walnut +\$62

91 white (CP) +\$62

Step 4. Pull Option/Finish


NN none +\$0

BK black +\$0

NH brushed nickel +\$0

## Storage Case Door

C5052.  
C5053.



### Product Information

#### Description

This door attaches to the front of a storage case. Attachment hardware is included.


#### Notes

Co/Struc® storage case drawers and storage case bottle drawers cannot be specified behind doors.

When specifying 17"- or 20"-high lockable door that has drawers above it, specify storage case shelf (C5041.) to prevent access to cabinet below.

For general information on locks, see Keyed-Alike Information in Appendices.

#### Dimensions



### Specification Information

#### Step 1.

**C505**

#### Step 2. Hinge Location

2. left-hand hinges
3. right-hand hinges

#### Step 3. Height

- 1722** 17" high  
**2022** 20" high  
**2522** 25" high

#### Step 4. Lock

- N** no lock  
**L** keyed differently  
**K** keyed alike

#### Prices for Steps 1-4.

		N	L	K
<b>C5052.</b>	<b>1722</b>	\$352	439	439
	<b>2022</b>	\$381	469	469
	<b>2522</b>	\$477	567	567
<b>C5053.</b>	<b>1722</b>	\$363	455	455
	<b>2022</b>	\$391	485	485
	<b>2522</b>	\$494	586	586


#### Step 5. Surface Finish

<b>91</b>	white (CP)	+\$31
<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

#### Step 6. Pull Option/Finish

<b>NN</b>	none	+\$0
<b>BK</b>	black	+\$0
<b>NH</b>	brushed nickel	+\$0



Storage Case LegC5050.




Product Information
Description
This leg supports a storage case mounted under a heavy-duty work surface. It provides support when the work surface holds heavy loads. The legs adjust in 1/2" increments to allow a work surface height of 34" to 38". Leg finish is medium tone.
Notes
Order storage case (C5040.) separately.
Dimensions

Specification Information
Step 1.
C5050.
\$278


Co/Struc® Storage








Product Information
Description
This unit hangs on a module, support panel, or wall strips above a work surface. It has sliding glass doors. Shipped knocked down.
Notes
Order optional task light separately: <ul style="list-style-type: none"><li>• Under shelf LED task light (G6170.)</li></ul>
Dimensions



Specification Information
Step 1.
C5014.16
Step 2. Width
30 30" wide
48 48" wide
Prices for Steps 1-2.
C5014.16 30 \$2010
48 \$2350
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
LU soft white +\$0
WN warm grey neutral +\$0
91 white (CP) +\$113



Shelf Unit

C5020.  
C5021.




Product Information
<div>Description</div> <p>This 3- or 5-shelf unit hangs from a module, panel, or wall strips. It has 9"- or 16"-deep shelves that adjust in 2½" increments. Shelves are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.</p>
<div>Notes</div> <p>For additional shelf, order add-on shelf (C5025.) separately.</p> <p>To cover shelf unit with dust cover, order add-on shelf (C5025.) and place in top position of unit.</p> <p>To enclose back of unit, order back panel (C5026.) separately.</p> <p>To hang unit on architectural wall and enclose back of unit, order shelf unit wall attachment panel (C5036.) separately.</p> <p>To enclose front of unit, order shelf unit door (C5027., C5028., C5032., C5037., or C5038.) separately.</p> <p>To ease cleaning and prevent the stacking of materials on top of shelf unit, order sloped shelf unit top (C0397.) separately.</p> <p>Order optional task light for 16"-deep unit separately:</p> <ul style="list-style-type: none"> <li>• Under shelf LED task light (G6170.)</li> </ul> <p>To attach task light to slanted shelf, order light bracket (C5051.) separately.</p>

Dimensions



3 Shelves



5 Shelves

**Specification Information****Step 1.****C502****Step 2. Shelves/Height**

- 0. 3 shelves and 24" high
- 1. 5 shelves and 34" high

**Step 3. Depth**

- 09 9" deep
- 16 16" deep

**Step 4. Width**

- 24 24" wide
- 30 30" wide
- 36 36" wide
- 48 48" wide

**Prices for Steps 1-4.**


	24	30	36	48
<b>C5020. 09</b>	\$698	738	753	793
<b>16</b>	\$874	920	971	1087
<b>C5021. 09</b>	\$980	1013	1055	1124
<b>16</b>	\$1176	1253	1329	1470

**Step 5. Side Panel/Shelf Finish**

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0
<b>91</b>	white (CP)	+\$62

Shelf Unit, Diagonal Corner

C5022.



Product Information

Description

This unit hangs from 24"-wide panels connected by an Action Office® 2-way 90° connector or wall strips used in a corner position. Shelves adjust in 2½" increments and are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.


Notes

For additional shelf in center section, order 24"-wide add-on shelf (C5025.1624). For additional shelf in side section, order 14"-wide add-on shelf (C5025.1614).

To cover back of center unit, order back panel (C5026.) separately.

To attach task light to slanted center shelf, order light bracket (C5051.) separately.

Dimensions




Specification Information



Step 1.		
C502		
Step 2. Shelves/Height		
2.1648	3 shelves and 24" high	
Prices for Steps 1-2.		
C5022.1648		\$2727
Step 3. Side Panel/Shelf Finish		
91	white (CP)	+\$164
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Shelf Unit, 90° Corner

C5024.




Product Information
Description
This 72"-high shelf unit hangs from 24"- or 48"-wide panels connected by 1 or 2 Action Office® 2-way 90° connectors or hangs from wall strips used in 1 or 2 90° corners. Shelves adjust in 2½" increments and are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.
Notes
For additional shelf, order add-on shelf (C5025.1631) separately.
To cover shelf unit with dust cover, order add-on shelf (C5025.1631) and place in top position of unit.
Dimensions



Specification Information
Step 1.
C5024.1648
\$4985
Step 2. Side Panel/Shelf Finish
HF inner tone light
LT light tone
LU soft white
WN warm grey neutral
91 white (CP)
+\$0
+\$0
+\$0
+\$0
+\$308

Co/Struc® Storage




Product Information

Description

This shelf fits inside a shelf unit, diagonal corner shelf unit, or 90° corner shelf unit. It is used as an additional shelf or dust cover.

Dimensions



Specification Information

Step 1.

C5025.

Step 2. Depth

09 9" deep

16 16" deep

Step 3. Width

For 9" deep (09)

24 22 3/4" wide

30 28 3/4" wide

36 34 3/4" wide

48 46 3/4" wide

For 16" deep (16)

14 13 1/2" wide

24 22 3/4" wide

30 28 3/4" wide

31 30 1/2" wide

36 34 3/4" wide


48 46 3/4" wide

Prices for Steps 1-3.

	14	24	30	31	36	48
C5025. 09	—	\$143	156	—	163	187
16	\$151	185	189	196	198	214

Step 4. Shelf Finish


91	white (CP)	+\$11
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Product Information
Description
This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.
Dimensions

Specification Information
Step 1.
CO39827
\$203

Co/Struc® Storage



Shelf Unit Back PanelC5026.



Product Information


Description

This panel attaches to a 3- or 5-shelf unit or diagonal corner shelf unit and closes the back. The shelf unit's top and bottom shelves must be flat.

Notes

Order shelf unit (C5020. or C5021.) or diagonal corner shelf unit (C5022.) separately.

Dimensions



Specification Information

Step 1.

C5026.

Step 2. Height

24	24" high
34	34" high

Step 3. Width

24	22¾" wide
30	28¾" wide
36	34¾" wide
48	46¾" wide

Prices for Steps 1-3.


		24	30	36	48
C5026.	24	\$279	302	324	386
	34	\$334	365	393	457

Step 4. Surface Finish

HF	inner tone light				+\$0
LT	light tone				+\$0
LU	soft white				+\$0
WN	warm grey neutral				+\$0




Light BracketC5051.



Product Information
Description
These brackets fasten a task light to the bottom of an overhead cabinet, a flat or slanted shelf in a 3- or 5-shelf unit, or a center shelf in a diagonal corner shelf unit. Package contains 1 pair. Finish is black.
Notes
For 16"-deep shelf, order under shelf LED task light (G6170.) separately.
Dimensions


Specification Information
Step 1.
C5051.
Step 2. Usage
16 for 16"-deep shelf
Prices for Steps 1-2.
C5051. 16\$177

Co/Struc® Storage



Shelf Unit Door

C5027.  
C5028.



Product Information

Description

This 24"-wide door encloses a shelf unit. Attachment hardware is included.

Notes

For 24"-high door, order 3-shelf unit (C5020.) separately.

For 34"-high door, order 5-shelf unit (C5021.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions

Specification Information

Step 1.

C502

Step 2. Hinge Location

- 7. left-hand hinges
- 8. right-hand hinges

Step 3. Height

- 2424 24" high
- 3424 34" high

Step 4. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-4.


		N	L	K
C5027.	2424	\$554	711	711
	3424	\$726	885	885
C5028.	2424	\$541	701	701
	3424	\$719	873	873

Step 5. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
91	white (CP)	+\$41


Step 6. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
NH	brushed nickel	+\$0



# Shelf Unit Double Doors

C5032.



## Product Information

### Description

These doors enclose a shelf unit. The doors are available with or without a lock. Attachment hardware is included.


### Notes

For 24"-high doors, order 3-shelf unit (C5020.) separately.

For 34"-high doors, order 5-shelf unit (C5021.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

### Dimensions



## Specification Information

### Step 1.

**C5032.**

### Step 2. Height

**24** 24" high

**34** 34" high

### Step 3. Width

**30** 30" wide

**36** 36" wide

**48** 48" wide

### Step 4. Lock

**N** no lock

**L** keyed differently

**K** keyed alike

### Prices for Steps 1-4.

			N	L	K
<b>C5032.</b>	<b>24</b>	<b>30</b>	\$598	767	767
		<b>36</b>	\$701	868	868
		<b>48</b>	\$859	1022	1022
	<b>34</b>	<b>30</b>	\$777	940	940
		<b>36</b>	\$912	1078	1078
		<b>48</b>	\$1193	1359	1359

### Step 5. Surface Finish


<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0
<b>91</b>	white (CP)	+\$52

### Step 6. Pull Option/Finish

<b>NN</b>	none	+\$0
<b>BK</b>	black	+\$0
<b>NH</b>	brushed nickel	+\$0

Sloped Shelf Unit Top

CO397.



Product Information

Description

This top attaches to a shelf unit. It is sloped to ease cleaning and prevent the stacking of materials on top of the cabinet.

Notes

Order shelf unit (C5020. or C5021.) separately.

Dimensions

Specification Information

Step 1.

CO397.

Step 2. Depth

09	10½" deep
16	17" deep


Step 3. Width

24	24" wide
30	30" wide
36	36" wide
48	48" wide

Prices for Steps 1-3.				
	24	30	36	48
CO397. 09	\$419	469	509	555
16	\$469	502	555	661


Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Shelf Unit Wall Attachment Panel

C5036.




Product Information
Description
This panel encloses the back of a 3- or 5-shelf unit and allows the shelf unit to hang on an architectural wall.
Notes
Order 3-shelf unit (C5020.) or 5-shelf unit (C5021.) separately.
Dimensions



Specification Information				
Step 1.				
C5036.				
Step 2. Height				
24	24" high			
34	34" high			
Step 3. Width				
24	24" wide			
30	30" wide			
36	36" wide			
48	48" wide			
Prices for Steps 1-3.				
	24	30	36	48
C5036. 24	\$334	372	450	522
34	\$361	387	457	545
Step 4. Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0
LU	soft white			+\$0
WN	warm grey neutral			+\$0
91	white (CP)			+\$26

B-Style Rail-Hanging Shelf

CO521.



Product Information

Description

This 15½"-high shelf with a lip hangs from a standard rail or adapter rail. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

For fully enclosed storage, order A/B-style flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately.


Order optional add-on shelves of equal depth and width separately:

- Add-on shelf with lip (C3019.)

Order optional task light separately:

- Under shelf LED task light (G6170.)

Dimensions



Specification Information

Step 1.

CO521.

Step 2. Depth

13	13" deep
16	16" deep

Step 3. Width

24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide

Prices for Steps 1-3.


	24	30	36	42	48
CO521. 13	\$615	631	656	673	710
16	\$849	868	912	933	1012

Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

## B-Style Shelf

C3011.



### Product Information

#### Description

This 15½"-high shelf with a lip hangs from a panel, a freestanding module, or wall strips. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

#### Notes

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

For fully enclosed storage of 24"- to 48"-wide shelf, order A-/B-style flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately.

Order optional add-on shelves of equal depth and width separately:


- Add-on shelf with lip (C3019.)

For 60"-wide shelf, order 60"-wide task light.

Order optional task light separately:

- Under shelf LED task light (G6170.)

### Dimensions



### Specification Information

#### Step 1.

**C3011.**

#### Step 2. Depth

<b>13</b>	13" deep
<b>16</b>	16" deep

#### Step 3. Width

<b>24</b>	24" wide
<b>30</b>	30" wide
<b>36</b>	36" wide
<b>42</b>	42" wide
<b>48</b>	48" wide
<b>60</b>	60" wide


#### Prices for Steps 1-3.

		<b>24</b>	<b>30</b>	<b>36</b>	<b>42</b>	<b>48</b>	<b>60</b>
<b>C3011. 13</b>		\$340	356	378	391	415	457
<b>16</b>		\$486	504	530	551	603	724

#### Step 4. Surface Finish

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

B-Style Shelf with Sloped TopC3091.



Product Information

Description

This 16"-deep shelf hangs from a panel, a freestanding module, or wall strips and has a sloped metal top. A filler piece is included for enclosed storage applications. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes



For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

To enclose back of shelf, order A-/B-style flipper door back panel (A3390.) separately.

Order optional task light separately:

- Under shelf LED task light (G6170.)

Dimensions



Specification Information

Step 1.

C3091.16

Step 2. Width

24	24" wide
30	30" wide
42	42" wide
48	48" wide

Prices for Steps 1-2.

C3091.16	24	\$674
	30	\$706
	42	\$821
	48	\$857


Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0



# Add-On Shelf with Lip

C3019.



## Product Information

### Description


This shelf fits into a B-style shelf or B-style rail-hanging shelf to add a second shelf. Attachment hardware is included.

### Notes

Order support components separately:

- B-style (C3011.) shelf
- B-style (CO521.) rail-hanging shelf

### Dimensions



## Specification Information

### Step 1.

**C3019.**

### Step 2. Depth

**13** 13" deep

**16** 16" deep

### Step 3. Width

**24** 24" wide

**30** 30" wide

**36** 36" wide

**42** 42" wide

**48** 48" wide

**60** 60" wide

### Prices for Steps 1-3.


		24	30	36	42	48	60
<b>C3019.</b>	<b>13</b>	\$198	208	231	245	266	298
	<b>16</b>	\$291	307	332	351	408	510

### Step 4. Surface Finish

<b>HF</b>	inner tone light						+\$0
<b>LT</b>	light tone						+\$0
<b>LU</b>	soft white						+\$0
<b>WN</b>	warm grey neutral						+\$0

Shelf Label Clip

CO398



Product Information

Description

This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.

Dimensions


Specification Information


Step 1.

CO39827

\$203

Co/Struc® Storage





### Product Information

#### Description

This enclosed storage unit hangs on rails for stationary storage. The locker has 20 drawer slots in 3" increments and holds interchangeable components. The locker has the following lock options: no lock, keyed lock, or keyless lock. The keyless lock mounts on the right-hand side and front of the locker. It includes a back-up keyed lock. The keyless lock is shipped with the battery.

#### Notes

Order rail components separately:


- Standard rail assembly (CO282)

Order transport/storage components separately:

- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)

For general information on locks, see Keyed-Alike Information in Appendices.

#### Dimensions



### Specification Information

#### Step 1.

##### CO561

#### Step 2. Lock

FF	no lock
FL	keyed lock, keyed differently
KA	keyed lock, keyed alike
EL	keyless lock, keyed differently
EK	keyless lock, keyed alike

#### Prices for Steps 1-2.

CO561	FF	\$2448
	FL	\$2667
	KA	\$2667
	EL	\$3943
	EK	\$3943

#### Step 3. Locker Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

#### Step 4. Door Finish

3B	beachglass	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Key Number


*For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK)*

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0

# Locker, with Tambour Door *continued*

Co/Struc® Storage

13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			



### Product Information

#### Description

This enclosed storage unit hangs on rails for stationary storage. The locker has 20 drawer slots in 3" increments and holds interchangeable components.

Door is hinged at the right.

Available with no lock, keyed lock, keyless lock, or RFID lock.

#### Notes

Door frame color matches locker color.

Order rail components separately:


- Standard rail assembly (CO282)

Order transport/storage components separately:

- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)

For general information on locks, see Keyed-Alike Information in Appendices.

#### Dimensions



### Specification Information

#### Step 1.

##### CO570

#### Step 2. Lock

FF	no lock
FL	keyed lock, keyed differently
KA	keyed lock, keyed alike
KP	keyless lock
RF	RFID lock

#### Prices for Steps 1-2.

CO570	FF	\$3160
	FL	\$3409
	KA	\$3409
	KP	\$3939
	RF	\$4900

#### Step 3. Locker Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

#### Step 4. Door Insert Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Key Number


##### For keyed lock, keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0

# Locker, with Hinged Door *continued*

Co/Struc® Storage

15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			




Product Information

**Description**  
This storage unit hangs on rails for stationary storage. The locker has 20 drawer slots in 3" increments and holds interchangeable components.

- Notes**
- Order rail components separately:
- Standard rail assembly (CO282)
- Order transport/storage components separately:
- C tray/shelf (CO203)
  - C wire shelf (CO252)
  - Cassette assembly (CO565)
  - Drawer (CO207, CO208, or CO209)
  - Pullout catheter rack (CO469.)
- Order optional locker tambour door (CO219) separately.

Dimensions



Specification Information


**Step 1.**  
**CO562FF** \$2161

**Step 2. Surface Finish**

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Locker on Wheel Base, with Tambour Door

CO564



Product Information

Description


This enclosed, mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The tambour door slides vertically for access to locker materials. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.

The locker has the following lock options: no lock, keyed lock, or keyless lock. The keyless lock mounts on the right-hand side and front of the locker. It includes a back-up keyed lock. The keyless lock is shipped with the battery.

Notes

- Order internal storage components separately:
- A-, B-, or C-size drawers (CO207, CO208, or CO209)
  - C tray/shelf (CO203)
  - C wire shelf (CO252)
  - Cassette assembly (CO565)
  - Pullout catheter rack (CO469.)
- For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions





Locker on Wheel Base, with Tambour  
Door *continued*

**Specification Information**

Step 1.

**CO564**

Step 2. Lock

<b>FF</b>	no lock
<b>FL</b>	keyed lock, keyed differently
<b>KA</b>	keyed lock, keyed alike
<b>EL</b>	keyless lock, keyed differently
<b>EK</b>	keyless lock, keyed alike

Prices for Steps 1-2.

<b>CO564 FF</b>	\$3417
<b>FL</b>	\$3690
<b>KA</b>	\$3585
<b>EL</b>	\$4911
<b>EK</b>	\$4911

Step 3. Locker Finish

<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>WN</b>	warm grey neutral	+\$0

Step 4. Door Finish

<b>3B</b>	beachglass	+\$0
<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>TR</b>	clear	+\$0
<b>WN</b>	warm grey neutral	+\$0

Step 5. Key Number

*For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK)*

<b>01</b>	key number 01	+\$0
<b>02</b>	key number 02	+\$0
<b>03</b>	key number 03	+\$0
<b>04</b>	key number 04	+\$0
<b>05</b>	key number 05	+\$0
<b>06</b>	key number 06	+\$0
<b>07</b>	key number 07	+\$0
<b>08</b>	key number 08	+\$0
<b>09</b>	key number 09	+\$0
<b>10</b>	key number 10	+\$0
<b>11</b>	key number 11	+\$0
<b>12</b>	key number 12	+\$0

<b>13</b>	key number 13	+\$0
<b>14</b>	key number 14	+\$0
<b>15</b>	key number 15	+\$0
<b>16</b>	key number 16	+\$0
<b>17</b>	key number 17	+\$0
<b>18</b>	key number 18	+\$0
<b>19</b>	key number 19	+\$0
<b>20</b>	key number 20	+\$0
<b>21</b>	key number 21	+\$0
<b>22</b>	key number 22	+\$0
<b>23</b>	key number 23	+\$0
<b>24</b>	key number 24	+\$0
<b>25</b>	key number 25	+\$0
<b>26</b>	key number 26	+\$0
<b>27</b>	key number 27	+\$0
<b>28</b>	key number 28	+\$0
<b>29</b>	key number 29	+\$0
<b>30</b>	key number 30	+\$0
<b>31</b>	key number 31	+\$0
<b>32</b>	key number 32	+\$0
<b>33</b>	key number 33	+\$0
<b>34</b>	key number 34	+\$0
<b>35</b>	key number 35	+\$0
<b>36</b>	key number 36	+\$0
<b>37</b>	key number 37	+\$0
<b>38</b>	key number 38	+\$0
<b>39</b>	key number 39	+\$0
<b>40</b>	key number 40	+\$0
<b>41</b>	key number 41	+\$0
<b>42</b>	key number 42	+\$0
<b>43</b>	key number 43	+\$0
<b>44</b>	key number 44	+\$0
<b>45</b>	key number 45	+\$0
<b>46</b>	key number 46	+\$0
<b>47</b>	key number 47	+\$0
<b>48</b>	key number 48	+\$0
<b>49</b>	key number 49	+\$0
<b>50</b>	key number 50	+\$0
<b>51</b>	key number 51	+\$0
<b>52</b>	key number 52	+\$0
<b>53</b>	key number 53	+\$0
<b>54</b>	key number 54	+\$0
<b>55</b>	key number 55	+\$0
<b>56</b>	key number 56	+\$0
<b>57</b>	key number 57	+\$0
<b>58</b>	key number 58	+\$0
<b>59</b>	key number 59	+\$0

Co/Struc® Storage


Locker on Wheel Base, with Tambour  
Door *continued*

Co/Struc® Storage

60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0
67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0
90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0

# Locker on Wheel Base, with Hinged Door

## CO574



### Product Information

#### Description

This enclosed, mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The door is hinged at the right and swings open for access to locker materials. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.

Available with no lock, keyed lock, keyless lock, or RFID lock.

#### Notes


Door frame color matches locker color.

Order internal storage components separately:

- A-, B-, or C-size drawers (CO207, CO208, or CO209)
- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Pullout catheter rack (CO469.)

For general information on locks, see Keyed-Alike Information in Appendices.

#### Dimensions



### Specification Information

#### Step 1.

#### CO574

#### Step 2. Lock

FF	no lock
FL	keyed lock, keyed differently
KA	keyed lock, keyed alike
KP	keyless lock
RF	RFID lock

#### Prices for Steps 1-2.

CO574	FF	\$4326
	FL	\$4572
	KA	\$4572
	KP	\$5104
	RF	\$6020

#### Step 3. Locker Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

#### Step 4. Door Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Key Number


#### For keyed lock, keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0

Locker on Wheel Base, with Hinged  
Door *continued*

Co/Struc® Storage

15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			




Product Information

**Description**  
This mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.

**Notes**  
Order internal storage components separately:  
• A-, B-, or C-size drawers (CO207, CO208, or CO209)  
• C tray/shelf (CO203)  
• C wire shelf (CO252)  
• Cassette assembly (CO565)  
• Pullout catheter rack (CO469.)  
Order optional locker tambour door (CO219) and locker lock kit (CO277) separately.  
Order optional locker hinged door (CO290) separately.


Dimensions



Specification Information

Step 1.		
CO563FF		\$3378
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Storage



Product Information

Description

This vertical door provides closure on a locker. It opens by sliding into the locker bottom.

Specification Information


Step 1.

CO219FF

\$372


Dimensions

Co/Struc® Storage



66-1/4" (1683)

SIDE




22-1/8" (562)

FRONT

Step 2. Surface Finish

3B	beachglass	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0




### Product Information

#### Description

This door provides closure on a locker. It is hinged at the right. Attachment hardware included.

Available with no lock, keyed lock, keyless lock or RFID lock.

#### Dimensions



### Specification Information

#### Step 1.

#### CO290

#### Step 2. Lock

FF	no lock
FL	keyed lock, keyed differently
KA	keyed lock, keyed alike
KP	keyless lock
RF	RFID lock

#### Prices for Steps 1-2.

CO290	FF	\$1004
	FL	\$1258
	KA	\$1258
	KP	\$1792
	RF	\$2805

#### Step 3. Door Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

#### Step 4. Insert Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Key Number

#### For keyed lock, keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0

# Locker Hinged Door *continued*

Co/Struc® Storage

15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			






# Locker Lock Kit *continued*

## Co/Struc® Storage

28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			
74	key number 74	+\$0			



Product Information


Description

This keyless lock converts a locker with tambour door to a lockable storage unit. The keyless lock mounts on the right-hand side and front of the locker. It includes a back-up keyed lock. The keyless lock is shipped with the battery.

Notes

For general information on locks, see Keyed-Alike Information in Appendices. Kit must be field installed.

Dimensions



Specification Information

Step 1.

CO279E

Step 2. Lock

K	keyless lock, keyed alike
L	keyless lock, keyed differently

Prices for Steps 1-2.

CO279E K	\$1601
L	\$1601

Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Step 4. Key Number


For keyless lock, keyed alike (K)

00	key number 00	+\$0
01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0

# Locker Lock Kit, Keyless *continued*

## Co/Struc® Storage


27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			




Product Information
Description
This component provides storage in an endoscopy/catheter cabinet. It has 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.
Notes
Tray/shelf surface with raised edges can be configured with subcontainers; order subcontainers separately.
Dimensions

Specification Information		
Step 1.		
CO204FF		\$121
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Storage





Co/Struc® Storage


Product Information

**Description**  
This frame hangs on a standard rail or adapter rail, or it is supported by an L cart, process table, counter top, or heavy-duty storage work surface. It has 4 slots at 3" intervals to enclose and support storage components and accessories. The C frame cannot mount under a 24"-wide work surface and it is not recommended for overhead storage.

Notes

- Order support components separately:
- Adapter rail (C1610.)
  - Counter top (CO555, or CO556)
  - L cart (CO212)
  - Process table (C4115., C4215., or C4313.)
  - Lab work table (C461G., C462G., or C462C.)
  - Standard rail (CO282)
  - Work table (C4312.)
  - Workstation (C4310. or C4311.)
- Order storage components separately:
- C tray/shelf (CO203)
  - Drawer (CO207, CO208, or CO209)

Dimensions



Specification Information


Step 1.	
CO205FF	\$321

Step 2. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Counter Top


CO555  
CO556




Product Information
Description
This surface mounts on a C frame to provide a work surface or storage surface. It includes a gripper package for attaching the top to a C frame. The single or extended counter top fits on 1 C frame. The extended counter top extends 4" past each C frame side.
Notes
Order C frame (CO205) separately.
Dimensions

Specification Information
Step 1.
CO
Step 2. Configuration/Width
555FF single
556FF extended
Prices for Steps 1-2.
CO555FF
CO556FF
Step 3. Surface Finish
HF inner tone light
LT light tone
LU soft white
WN warm grey neutral


Co/Struc® Storage




TOP



SINGLE  
FRONT



EXTENDED  
FRONT



Product Information

Description

This lid provides a top for a C frame. It has raised side edges to keep items from sliding off.

Notes

Use lid when storing L cart below 38"-high process table or cantilevered work surface; existing drawer bearers must be removed.

Dimensions

Specification Information


Step 1.

CO201FF	\$225
---------	-------


Step 2. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Storage







**Product Information**

**Description**

This component provides storage within a C frame, locker, or interior cabinet adapter. It is designed with 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.

**Notes**

Order support component separately:


- C frame (CO205)
- Interior cabinet adapter (CG591.)
- Locker or locker on wheel base (CO561, CO562, CO563, or CO564)

Tray/shelf surface with raised edges can be configured with subcontainers or subdivider vane. Order subdivider vane (CO242) or subcontainers separately.


**Dimensions**

Specification Information		
Step 1.		
CO203FF		\$152
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0


Co/Struc® Storage



TOP



FRONT



Product Information

Description

This component organizes IV bags, forms, tape rolls, and other small items within a drawer or C tray/shelf.

Dimensions

Specification Information


Step 1.


**CO242FF** \$75

Step 2. Surface Finish

**LU** soft white +\$0

Co/Struc® Storage







Product Information
<b>Description</b>
This wire shelf fits into the slots of a C frame, locker, or interior cabinet adapter and stores items requiring air circulation. It cannot be used with a subdivider vane or subdivider.
<b>Notes</b>
Order support component separately:
• C frame (CO205)
• Interior cabinet adapter (CG591.)
• Locker or locker on wheel base (CO561, CO562, CO563, or CO564)
<b>Dimensions</b>

Specification Information		
Step 1.		
CO252FF		\$195
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage






Product Information	
Description	This hardware package mounts a C frame under a work surface. 1 gripper package is required to mount each C frame to a work surface. Attachment screws are included.
Notes	C frame cannot mount with grippers under heavy-duty work surface.
Dimensions	


Specification Information		
Step 1.		
CO266FF		\$94
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Co/Struc® Storage



## Drawer

CO207  
CO208  
CO209



### Product Information

#### Description

This drawer fits into a locker or C frame and has stops to prevent accidental removal. It can also be mounted under a work surface using drawer bearers, ordered separately. The drawer is available with no lock, keyed lock, keyless lock or RFID lock.

#### Notes

To mount drawer under table or work surface, order drawer bearer package (CO228) or lockable drawer bearer (CO396) separately.


Order optional A-size subcontainer (CO231, CO232, CO233, CO234, CO235, or CO369) or B-size subcontainer (CO236, CO238, or CO240) separately.

Drawer with keyless lock (KP) or RFID (RF) option cannot be used in locker with tambour door.


Keyed lock is located on right side; keyless lock and RFID lock are located on left side.

For general information on locks, see Keyed-Alike Information in Appendices.


### Dimensions




TOP



A SIZE  
FRONT



B SIZE  
FRONT



C SIZE  
FRONT

**Specification Information**

Step 1.

**CO20**

Step 2. Size/Height

<b>7</b>	a size	
<b>8</b>	b size	
<b>9</b>	c size	

Step 3. Lock

<b>FF</b>	no lock	
<b>FL</b>	keyed differently	
<b>KA</b>	keyed alike	
<b>KP</b>	keyless lock	

Prices for Steps 1-3.

<b>CO207</b>	<b>FF</b>	\$139
	<b>FL</b>	\$290
	<b>KA</b>	\$290
	<b>KP</b>	\$924

<b>CO208</b>	<b>FF</b>	\$180
	<b>FL</b>	\$328
	<b>KA</b>	\$328
	<b>KP</b>	\$959

<b>CO209</b>	<b>FF</b>	\$225
	<b>FL</b>	\$376
	<b>KA</b>	\$376
	<b>KP</b>	\$1007

Step 4. Surface Finish

<b>2B</b>	juniper	+\$0
<b>3B</b>	beachglass	+\$0
<b>6F</b>	midnight	+\$0
<b>DA</b>	deep blue	+\$0
<b>HF</b>	inner tone light	+\$0
<b>LT</b>	light tone	+\$0
<b>LU</b>	soft white	+\$0
<b>VR</b>	vivid red	+\$0
<b>WN</b>	warm grey neutral	+\$0
<b>YT</b>	yellow	+\$0

Step 5. Key Number

*For keyed alike (KA)*

<b>01</b>	key number 01	+\$0
<b>02</b>	key number 02	+\$0

<b>03</b>	key number 03	+\$0
<b>04</b>	key number 04	+\$0
<b>05</b>	key number 05	+\$0
<b>06</b>	key number 06	+\$0
<b>07</b>	key number 07	+\$0
<b>08</b>	key number 08	+\$0
<b>09</b>	key number 09	+\$0
<b>10</b>	key number 10	+\$0
<b>11</b>	key number 11	+\$0
<b>12</b>	key number 12	+\$0
<b>13</b>	key number 13	+\$0
<b>14</b>	key number 14	+\$0
<b>15</b>	key number 15	+\$0
<b>16</b>	key number 16	+\$0
<b>17</b>	key number 17	+\$0
<b>18</b>	key number 18	+\$0
<b>19</b>	key number 19	+\$0
<b>20</b>	key number 20	+\$0
<b>21</b>	key number 21	+\$0
<b>22</b>	key number 22	+\$0
<b>23</b>	key number 23	+\$0
<b>24</b>	key number 24	+\$0
<b>25</b>	key number 25	+\$0
<b>26</b>	key number 26	+\$0
<b>27</b>	key number 27	+\$0
<b>28</b>	key number 28	+\$0
<b>29</b>	key number 29	+\$0
<b>30</b>	key number 30	+\$0
<b>31</b>	key number 31	+\$0
<b>32</b>	key number 32	+\$0
<b>33</b>	key number 33	+\$0
<b>34</b>	key number 34	+\$0
<b>35</b>	key number 35	+\$0
<b>36</b>	key number 36	+\$0
<b>37</b>	key number 37	+\$0
<b>38</b>	key number 38	+\$0
<b>39</b>	key number 39	+\$0
<b>40</b>	key number 40	+\$0
<b>41</b>	key number 41	+\$0
<b>42</b>	key number 42	+\$0
<b>43</b>	key number 43	+\$0
<b>44</b>	key number 44	+\$0
<b>45</b>	key number 45	+\$0
<b>46</b>	key number 46	+\$0
<b>47</b>	key number 47	+\$0
<b>48</b>	key number 48	+\$0
<b>49</b>	key number 49	+\$0

# Drawer *continued*

50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			
74	key number 74	+\$0			
75	key number 75	+\$0			
76	key number 76	+\$0			
77	key number 77	+\$0			
78	key number 78	+\$0			
79	key number 79	+\$0			
80	key number 80	+\$0			
81	key number 81	+\$0			
82	key number 82	+\$0			
83	key number 83	+\$0			
84	key number 84	+\$0			
85	key number 85	+\$0			
86	key number 86	+\$0			
87	key number 87	+\$0			
88	key number 88	+\$0			
89	key number 89	+\$0			
90	key number 90	+\$0			
91	key number 91	+\$0			
92	key number 92	+\$0			
93	key number 93	+\$0			
94	key number 94	+\$0			
95	key number 95	+\$0			
96	key number 96	+\$0			


Co/Struc® Storage

## Drawer with Lockable Lid

CO370

CO372

CO374



### Product Information

#### Description

This drawer has a hinged steel lid that can be locked. The lid can open when the drawer is pulled out to a stop position; the lid cannot be removed.

The drawer is 3", 6", or 9" high.


#### Notes

For B-size drawer with lockable lid, order A-size subcontainer only. For


C-size drawer with lockable lid, order A- or B-size subcontainers only.

For general information on locks, see Keyed-Alike Information in Appendices.


#### Dimensions




TOP



A SIZE  
FRONT



B SIZE  
FRONT



C SIZE  
FRONT

### Specification Information

#### Step 1.

#### CO37

#### Step 2. Size/Height

0	a size
2	b size
4	c size

#### Step 3. Lid Lock

FF	keyed differently
KA	keyed alike

#### Prices for Steps 1-3.

CO370	FF	\$491
	KA	\$491
CO372	FF	\$529
	KA	\$529
CO374	FF	\$581
	KA	\$581

#### Step 4. Drawer Finish

2B	juniper	+\$0
3B	beachglass	+\$0
6F	midnight	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
VR	vivid red	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Lid Finish

LU	soft white	+\$0
----	------------	------

#### Step 6. Key Number

##### For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0




# Drawer with Lockable Lid *continued*

09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			

Co/Struc® Storage

## Locked Drawer with Lockable Lid

CO371  
CO373  
CO375



### Product Information

#### Description

This drawer fits into a locker and has stops to prevent accidental removal. It has a lockable, non-removable lid that can open when the drawer is pulled out. The drawer is available with a keyed lock, keyless lock or RFID lock.

#### Notes


Locked drawer with lockable lid cannot be used with a C frame.

Drawer with keyless lock (KP) or RFID lock (RF) option cannot be used in locker with tambour door.




Keyed lock is located on right side; keyless lock and RFID lock are located on left side.

For general information on locks, see Keyed-Alike Information in Appendices.

### Dimensions



TOP



### Specification Information

#### Step 1.

#### CO37

#### Step 2. Size/Height

1	a size
3	b size
5	c size

#### Step 3. Drawer/Lid Lock

FF	keyed differently
KA	keyed alike
KP	keyless lock

#### Prices for Steps 1-3.

CO371	FF	\$632
	KA	\$632
	KP	\$1290
CO373	FF	\$672
	KA	\$672
	KP	\$1336
CO375	FF	\$722
	KA	\$740
	KP	\$1374

#### Step 4. Drawer Finish

2B	juniper	+\$0
3B	beachglass	+\$0
6F	midnight	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
VR	vivid red	+\$0
WN	warm grey neutral	+\$0

#### Step 5. Drawer Key Number

##### For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0

# Locked Drawer with Lockable Lid

continued

08	key number 08	+\$0	55	key number 55	+\$0
09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			

Co/Struc® Storage

## Locked Drawer with Lockable Lid

continued

Co/Struc® Storage

Step 6. Lid Finish		
LU	soft white	+\$0

Step 7. Lid Key Number		
<i>For keyed alike (KA) or keyless lock (KP)</i>		


01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0

43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0
50	key number 50	+\$0
51	key number 51	+\$0
52	key number 52	+\$0
53	key number 53	+\$0
54	key number 54	+\$0
55	key number 55	+\$0
56	key number 56	+\$0
57	key number 57	+\$0
58	key number 58	+\$0
59	key number 59	+\$0
60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0
67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0

## Locked Drawer with Lockable Lid

*continued*

90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0



Product Information

Description

This clear lid fits on a Co/Struc drawer and provides top closure.

Notes

Order drawer (CO207, CO208, or CO209) separately.

Dimensions


Specification Information


Step 1.

CO508.20

\$207

Co/Struc® Storage







Product Information
Description
This package is used to mount a nonlocking drawer under a work surface. The drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.
Dimensions

Specification Information		
Step 1.		
CO228FF		\$51
Step 2. Surface Finish		
BK	black	+\$0

Co/Struc® Storage






Product Information
Description
This component is used to mount a drawer with lock under a work surface. The lockable drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.
Dimensions

Specification Information		
Step 1.		
CO396FF		\$113
Step 2. Surface Finish		
BK	black	+\$0


Co/Struc® Storage






Label Clip

CO274  
CO462




Product Information
Description
This clip snaps onto a drawer, tray/shelf, or flipper to hold content labels.
The clip accepts pressure-sensitive labels from the front or slide-in labels from the back. Finish is clear. Package contains 12 clips; labels not included.
Dimensions



Specification Information
Step 1.
CO
Step 2. Width
27427 6" wide
462 19" wide
Prices for Steps 1-2.
CO27427 \$50
CO462 \$88



Co/Struc® Storage



Product Information
Description
This assembly fits inside a locker or C frame or fits below a work surface with a drawer bearer. It has a wire cassette, cassette bins, and subdividers. The assembly fills 2 slots in a locker or C frame.
Dimensions


Specification Information		
Step 1.		
CO565FF		\$748
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage



## Dispensing Rail

CO365  
CO366  
CO367  
CO368



### Product Information

#### Description

This rail attaches to a panel, a freestanding module, or wall strips to support subcontainers at a 15° angle for gravity-feed dispensing. Package contains 2.

The A-size rail holds A-size subcontainers; the B-size rail holds B-size subcontainers.

#### Notes

Dispensing rails hold the following number of subcontainers:

A-Size Rail Width—Number of Subcontainers

24"—5 A-1, A-2, A-3 or A-4

24"—2 A-5

24"—10 A-8

48"—10 A-1, A-2, A-3, or A-4

48"—5 A-5

48"—20 A-8

B-Size Rail Width—Number of Subcontainers

24"—5 B-1 or B-3


24"—2 B-5

48"—10 B-1 or B-3

48"—5 B-5

Order optional A-size subcontainer (CO231, CO232, CO233, CO234, CO235, or CO369) or B-size subcontainer (CO236, CO238, or CO240) separately.

### Dimensions



Co/Struc® Storage

Specification Information

Step 1.

CO36

Step 2. Size/Width

- 5FF a size and 24" wide
- 6FF a size and 48" wide
- 7FF b size and 24" wide
- 8FF b size and 48" wide


Prices for Steps 1-2.

CO365FF	\$520
CO366FF	\$664
CO367FF	\$560
CO368FF	\$668

Step 3. Surface Finish

LU soft white	+\$0
---------------	------

## Subcontainer



CO231  
CO232  
CO233  
CO234  
CO235  
CO236  
CO238  
CO240  
CO369

### Product Information

#### Description

This storage container fits into a drawer or onto a dispensing rail to organize and store small items. Package contains 12.

The A-size subcontainer is used with an A-size drawer (without lid) or A-size dispensing rail. It is  $2\frac{3}{16}$ " high; depth and width dimensions are:

Type—Depth—Width

A-1— $4\frac{9}{16}$ "— $4\frac{9}{16}$ "

A-2— $9\frac{3}{16}$ "— $4\frac{9}{16}$ "

A-3— $13\frac{13}{16}$ "— $4\frac{9}{16}$ "

A-4— $18\frac{7}{16}$ "— $4\frac{9}{16}$ "

A-5— $9\frac{3}{16}$ "— $9\frac{3}{16}$ "

A-8— $13\frac{13}{16}$ "— $2\frac{1}{4}$ "

The B-size subcontainer is used with a B-size drawer (without lid), B-size dispensing rail, or C-size drawer. It is  $5\frac{3}{16}$ " high; depth and width dimensions are:

Type—Depth—Width

B-1— $4\frac{9}{16}$ "— $4\frac{9}{16}$ "

B-3— $13\frac{13}{16}$ "— $4\frac{9}{16}$ "


B-5— $9\frac{3}{16}$ "— $9\frac{3}{16}$ "

#### Notes


Order following products separately:

- Drawer (CO207, CO208, or CO209)
- Dispensing rail (CO365, CO366, CO367, or CO368)

### Dimensions



TOP



A-SIZE  
FRONT

B-SIZE  
FRONT

Co/Struc® Storage

Specification Information

Step 1.

CO

Step 2. Size

- 231FF a-1
- 232FF a-2
- 233FF a-3
- 234FF a-4
- 235FF a-5
- 369FF a-8
- 236FF b-1
- 238FF b-3
- 240FF b-5

Prices for Steps 1-2.


CO231FF	\$113
CO232FF	\$113
CO233FF	\$123
CO234FF	\$123
CO235FF	\$123
CO369FF	\$118
CO236FF	\$123
CO238FF	\$169
CO240FF	\$169

Step 3. Surface Finish

LU	soft white	+\$0
----	------------	------

## A-Size Drawer Organizer Kit

CO540  
CO541  
CO542  
CO543  
CO544  
CO545



### Product Information

#### Description

This kit fits inside an A-size drawer (without lid) and has 6 combinations of subcontainers and subcontainer dividers.

A-size subcontainer is 2 $\frac{3}{16}$ " high. Width and depth dimensions are:

Type—Width—Depth

A-1—4 $\frac{9}{16}$ "—4 $\frac{9}{16}$ "

A-2—4 $\frac{9}{16}$ "—9 $\frac{3}{16}$ "



A-3—4 $\frac{9}{16}$ "—13 $\frac{13}{16}$ "

A-4—4 $\frac{9}{16}$ "—18 $\frac{7}{16}$ "

A-5—9 $\frac{3}{16}$ "—9 $\frac{3}{16}$ "


A-8—2 $\frac{1}{4}$ "—13 $\frac{13}{16}$ "

#### Dimensions




CO540FF


TOP




CO541FF




CO542FF



CO543FF



CO544FF



CO545FF

### Specification Information

#### Step 1.

#### CO54

#### Step 2. Components

0FF 4 a-3 subcontainers, 6 a-3 dividers

1FF 1 a-3 subcontainer, 6 a-8 subcontainers, 2 a-3 dividers, 7 a-8 dividers

2FF 2 a-1 subcontainers, 2 a-2 subcontainers, 1 a-3 subcontainer, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers

3FF 3 a-3 subcontainers, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers

4FF 4 a-2 subcontainers, 1 a-4 subcontainer, 8 a-3 dividers

5FF 2 a-1 subcontainers, 2 a-3 subcontainers, 1 a-5 subcontainer, 4 a-3 dividers

#### Prices for Steps 1-2.


CO540FF	\$104
CO541FF	\$110
CO542FF	\$110
CO543FF	\$104
CO544FF	\$104
CO545FF	\$104

#### Step 3. Surface Finish

LU	soft white	+\$0
----	------------	------

B-Size Drawer Organizer Kit



CO546  
CO547



Product Information

**Description**  
This kit fits inside a B-size drawer (without lid) or C-size drawer and has 2 combinations of subcontainers and subcontainer dividers.  
The B-size subcontainer is 5<sup>3</sup>/<sub>16</sub>" high. Width and depth dimensions are:  
Type—Width—Depth  
B-1—4<sup>9</sup>/<sub>16</sub>"—4<sup>9</sup>/<sub>16</sub>"  
B-3—4<sup>9</sup>/<sub>16</sub>"—13<sup>13</sup>/<sub>16</sub>"  
B-5—9<sup>3</sup>/<sub>16</sub>"—9<sup>3</sup>/<sub>16</sub>"

Dimensions




Specification Information

Step 1.		
CO54		
Step 2. Components		
6FF	4 b-3 subcontainers, 4 b-3 dividers	
7FF	2 b-1 subcontainers, 2 b-3 subcontainers, 1 b-5 subcontainer, 2 b-3 dividers	
Prices for Steps 1-2.		
CO546FF		\$125
CO547FF		\$125
Step 3. Surface Finish		
LU	soft white	+\$0




Subcontainer Divider

CO378  
CO379  
CO380



Product Information
Description
This divider fits into subcontainers to separate and organize stored items. A divider can be adjusted by sliding it to the desired position. The divider is not available for A-5 and B-5 subcontainers. Package contains 12.
The divider for an A-8 subcontainer is 1⅝" high and 2⅞" wide; the divider for an A-1, A-2, A-3, or A-4 subcontainer is 1⅝" high and 4¼" wide; and the divider for a B-1 or B-3 subcontainer is 4" high and 4¼" wide.
Dimensions




Specification Information
Step 1.
CO3
Step 2. Usage
78FF for a-8 subcontainer
79FF for a-1, a-2, a-3, or a-4 subcontainers
80FF for b-1 or b-3 subcontainers
Prices for Steps 1-2.
CO378FF \$106
CO379FF \$106
CO380FF \$108
Step 3. Surface Finish
LU soft white +\$0

Co/Struc® Storage

Subcontainer Handle

C0376  
C0377



Product Information

Description

This snap-on handle has a surface for identification labels in a picking station. It also can convert a small or large subcontainer into a drawer when used on the flat surface of a C tray. The handle cannot be used with a subcontainer in an A-size drawer. Package contains 12.

Notes

Handle accepts following label:

- 2 1/8"-wide handle: 3/4" x 1" gummed label
- 4 1/2"-wide handle: 3/4" x 3" drop-in label

For 4 1/2"-wide handle, order label (C0459) separately.

Dimensions

Specification Information

Step 1.

C037

Step 2. Width

6FF	2 1/8" wide
7FF	4 1/2" wide


Prices for Steps 1-2.


C0376FF	\$56
C0377FF	\$70

Step 3. Surface Finish

LU	soft white	+\$0
----	------------	------

Co/Struc® Storage





Product Information


**Description**  
This label slides into a 4½"-wide subcontainer handle. Package contains 100 cards with 10 perforated labels per card (1000 labels total).

Specification Information


**Step 1.**  
**CO459** \$47

Subcontainer Label Flag

C0381  
C0382



Product Information
Description
This flag snaps onto subcontainers that are hung from a dispensing rail. It provides a visible surface for identification labels. Package contains 12.
The flag for an A-8 subcontainer is 1¼" wide; the flag for an A-3 or B-3 subcontainer is 3½" wide.
Dimensions




Specification Information
Step 1.
C038
Step 2. Usage
1FF for a-8 subcontainer
2FF for a-3 or b-3 subcontainer
Prices for Steps 1-2.
C0381FF
C0382FF
Step 3. Surface Finish
LU soft white

\$152  
\$172


+\$0

Pullout Catheter Rack

CO469.




Product Information
Description
This 20"-wide, powder-coated wire rack fits into the slots of a locker to organize catheters. It is available with 3 or 4 independent, pullout slides. Each pullout slide has 4 prongs, 4 side label clips, and 1 front label holder. Each prong holds up to 10 catheters. Finish is soft white.
Dimensions



Specification Information
Step 1.
CO469.20
Step 2. Slides
03 3 slides
04 4 slides
Prices for Steps 1-2.
CO469.20 03 \$1091
04 \$1219

Co/Struc® Storage




Product Information

Description

This height-adjustable laptop cart delivers point-of-care technology and holds a 17" laptop. The surface provides security by locking a laptop into place during use. The cart holds a maximum of 30 pounds and has 4"-high casters; 2 locking and 2 nonlocking. It has an 18" height-adjustment range and is adjusted with a paddle on the front of the surface. The cart includes a swivel mouse tray. Shipped fully assembled.

The surface is Cloud Surf laminate; the base finish is polished aluminum.

Dimensions



Specification Information

Step 1.

CO600.

Step 2. Power Option

- N no power strip
- P power strip

Step 3. Bar Code Scanner Holder

- N no bar code scanner holder
- B bar code scanner holder

Step 4. Utility Basket

- N no utility basket
- U utility basket

Step 5. Drawer


- N no drawer
- D drawer

Prices for Steps 1-5.

				NN	ND	UN	UD
CO600.	N	N		\$3963	4531	4263	4830
		B		\$4116	4679	4416	4978
	P	N		\$4468	5030	4770	5332
		B		\$4616	5180	4917	5480

Flat Panel Cart

CO601.




**Product Information**

**Description**

This height-adjustable cart delivers point-of-care technology. It holds a flat panel monitor arm and includes a 10"-wide column-mounted CPU holder with straps. The arm meets VESA standards (Video Electronic Standards Association - Mounting Interface Standards). The cart holds a maximum of 30 pounds and has 4"-high casters; 2 locking and 2 nonlocking. It has an 18" height-adjustment range and is adjusted with a paddle on the front of the surface. The cart includes a swivel mouse tray and 2 plates: 75mm and 100mm. Shipped fully assembled.

The surface is Cloud Surf laminate; the base finish is polished aluminum; and the CPU holder is black.

**Dimensions**



**Specification Information**

**Step 1.**

**CO601.**

**Step 2. Power Option**

**N** no power strip

**P** power strip

**Step 3. Bar Code Scanner Holder**

**N** no bar code scanner holder

**B** bar code scanner holder

**Step 4. Utility Basket**

**N** no utility basket

**U** utility basket

**Step 5. Drawer**

**N** no drawer

**D** drawer


**Prices for Steps 1-5.**

			NN	ND	UN	UD
CO601.	N	N	\$4181	4777	4497	5094
		B	\$4341	4937	4660	5252
	P	N	\$4710	5303	5029	5623
		B	\$4871	5464	5188	5780

Co/Struc® Transport

Coiled Power Cord

CO602.



Product Information

Description

This hospital-grade power cord has a 1.5' coiled cord that extends to 7.5'. Finish is black.

Dimensions

Specification Information

Step 1.

CO602.

\$267


Co/Struc® Transport






Utility Basket


CO605.



Product Information
Description
This perforated basket attaches to the column of a laptop or flat-panel cart and holds a utility basket. It can be used to hold charts, paper products, and other supplies. Finish is black.
Dimensions

Specification Information
Step 1.
CO605.
\$338







Product Information
Description
This holder attaches to the surface of a laptop or flat-panel cart. It manages cables and wires and includes reusable cable ties. Finish is black.
Notes
Laptop-cart cable manager has 1 solid holder. Flat-panel-cart cable manager has 2 pieces that attach to each side of the flat-screen arm holder and go down through the work surface.
Dimensions

Specification Information		
Step 1.		
CO608.		
Step 2. Usage		
01	for laptop cart	
02	for flat panel cart	
Prices for Steps 1-2.		
CO608. 01		\$244
02		\$244


Co/Struc® Transport



01




02



Product Information
Description
This mobile cart has a platform, 4 rail profiles for hanging frames, and 4 casters (2 swivel and 2 swivel with brake) that are 4¾" hard-rubber zinc.
Notes
To hang drawer or shelf storage components, order C frame (CO205) separately.
Dimensions

Specification Information
Step 1.
CO212
Step 2. Application
FF standard
Prices for Steps 1-2.
CO212 FF
\$1138
Step 3. Surface Finish
LU soft white
+\$0

Co/Struc® Transport





## Index by Product Name

Co/Struc® System	
2-Way 90° Corner Connector	page(s) 31
2-Way Module Corner Connector	44
3-Way 90° Corner Connector	32
3-Way Module Corner Connector	45
Access Panel Assembly	52
Adapter Rail	7
ADA Sink Support, Systems Mounted	108
Add-On Shelf, Shelf Unit	126
Add-On Shelf with Lip	137
Adjustable Work Table Base	95
Anchor Support Kit	47
A-Size Drawer Organizer Kit	183
Backsplash Filler	79
Basic Work Table Base	90
B-Size Drawer Organizer Kit	184
B-Style Rail-Hanging Shelf	134
B-Style Shelf	135
B-Style Shelf with Sloped Top	136
Cable Manager	194
Cassette Assembly	178
Ceiling Entry Support Panel	26
Center Closure Panel	59
C Frame	158
C Frame Support	76
Clear Lid	174
C Lid	160
Coiled Power Cord	192
Combination Access Panel Assembly	54
Component Adapter Kit	77
Counter Top	159
CS Tray/Shelf	157
C Tray/Shelf	161
C Wire Shelf	163
Dispensing Rail	179
Drawer	165
Drawer Bearer Package	175
Drawer with Lockable Lid	168
Draw Rod, Seismic	12
End Cap for Standard or Seismic Rail	10
End Closure Assembly	58
Equipment Rail, Systems Mounted	9
Equipment Rail, Wall Mounted	8
Flat Panel Cart	191
Flat Panel Mount	105
Floor Anchor Adapter Kit	17
Floor Anchor Bracket, Action Office® Series 1	14

Floor Anchor Bracket, Action Office® Series 2	15
Floor Anchor Bracket, Co/Struc® Lab Module	13
Floor Entry Support Panel	24
Gripper Package	164
Heavy-Duty Corner Work Surface	80
Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2 90° Corner Connector	82
Heavy-Duty Work Surface	70
Interface Connector	33
Interior Filler	112
Interior Shelf	49
Interior Shelf and Support Kit	51
Label Clip	177
Laptop Cart	190
L Cart	195
Light Bracket	129
Lockable Drawer Bearer	113, 176
Locked Drawer with Lockable Lid	170
Locker	143
Locker Hinged Door	151
Locker Lock Kit	153
Locker Lock Kit, Keyless	155
Locker on Wheel Base	149
Locker on Wheel Base, with Hinged Door	147
Locker on Wheel Base, with Tambour Door	144
Locker Tambour Door	150
Locker, with Hinged Door	141
Locker, with Tambour Door	139
Low Cup Sink Access Panel Assembly	56
Low Module	38
Low Module Terminal Panel	41
Marker Board	104
Mid-Run Support, Double Sided	36
Mid-Run Support, Single Sided	35
Module Support Leg	46
Module Wall Start	43
Open Support Panel	20
Panel Connector	34
Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	18
Panel Wall Start	30
Pass-Through Cover	75
Pass-Through Module	39
Power Strip	103
Process Table, Electric Height Adjustment	83
Pullout Catheter Rack	189
Rail-Hanging Work Surface	74
Rail Hardware	6

# Index by Product Name *continued*

Index: Product Name

Reagent Shelf	84
Riser Frame Divider	101
Riser Frame End Cover	102
Seismic Finished End	16
Shelf End Cap	50
Shelf Label Clip	127, 138
Shelf Unit	122
Shelf Unit, 90° Corner	125
Shelf Unit Back Panel	128
Shelf Unit, Diagonal Corner	124
Shelf Unit Door	130
Shelf Unit Double Doors	131
Shelf Unit Wall Attachment Panel	133
Sink Module	106
Sink Unit for Storage	107
Sloped Shelf Unit Top	132
Standard-Duty Work Surface	68
Standard Rail	5
Standard Work Table Base	92
Storage Case	114
Storage Case Bottle Drawer	118
Storage Case Door	119
Storage Case Drawer	117
Storage Case Leg	120
Storage Case Pullout Shelf	116
Storage Case Shelf	115
Storage Unit	109
Storage Unit Door	111
Storage Unit Shelf	110
Storage Unit with Sliding Doors	121
Subcontainer	181
Subcontainer Divider	185
Subcontainer Handle	186
Subcontainer Handle Label	187
Subcontainer Label Flag	188
Subdivider Vane	162
Support Panel	22
Support Panel Top Cap	37
Terminal Panel	28
Top Shelf	48
Utility Basket	193
Vertical Chase Cover	42
Wall Rail, Seismic	11
Wall Strip Cover	19
Wire Chase Overlay	67
Wire Chase Rail	60
Wire Chase Rail, 90° Inside Corner	62

Wire Chase Rail Bridge	64
Wire Chase Rail Cover	66
Wire Chase Rail End Cap	65
Workstation, Electric Height Adjustment	85
Workstation, Manual Height Adjustment	87
Work Surface Backsplash	78
Work Surface Support	73
Work Table, Electric Height Adjustment	89
Work Table Reagent Shelf	100
Work Table Surface, Rectangular	97

## Index by Product Number

A1250. Seismic Finished End	page(s) 16
C460G. Basic Work Table Base	90
C461G. Standard Work Table Base	92
C462C. Adjustable Work Table Base	95
C462G.	
C1112. Low Module	38
C1114. Pass-Through Module	39
C1122. Low Module Terminal Panel	41
C1150. Open Support Panel	20
C1151. Support Panel	22
C1155. Floor Entry Support Panel	24
C1156. Ceiling Entry Support Panel	26
C1161. Terminal Panel	28
C1170. Equipment Rail, Wall Mounted	8
C1172. Equipment Rail, Systems Mounted	9
C1211. Panel Wall Start	30
C1212. Module Wall Start	43
C1220. 2-Way Module Corner Connector	44
C1221. 2-Way 90° Corner Connector	31
C1230. 3-Way Module Corner Connector	45
C1231. 3-Way 90° Corner Connector	32
C1261. Interface Connector	33
C1272. Mid-Run Support, Single Sided	35
C1273. Mid-Run Support, Double Sided	36
C1274. Support Panel Top Cap	37
C1275. Module Support Leg	46
C1276. Anchor Support Kit	47
C1281. Panel Connector	34
C1313. Wire Chase Rail	60
C1314.	
C1315. Wire Chase Rail Cover	66
C1316.	
C1317.	
C1318.	
C1321. Wire Chase Rail, 90° Inside Corner	62
C1322.	
C1326. Wire Chase Rail Bridge	64
C1327.	
C1331. Wire Chase Rail End Cap	65
C1332.	
C1358. Wire Chase Overlay	67
C1412. Access Panel Assembly	52
C1422. Combination Access Panel Assembly	54
C1432. Low Cup Sink Access Panel Assembly	56
C1451. End Closure Assembly	58
C1455. Center Closure Panel	59
C1460. Vertical Chase Cover	42

C1510. Shelf End Cap	50
C1510. Top Shelf	48
C1511. Interior Shelf	49
C1511. Shelf End Cap	50
C1515. Interior Shelf and Support Kit	51
C1610. Adapter Rail	7
C2140. Work Surface Support	73
C2141. Heavy-Duty Work Surface	70
C2142. Standard-Duty Work Surface	68
C2192. C Frame Support	76
C2193. Component Adapter Kit	77
C2215. Heavy-Duty Corner Work Surface	80
C2216. Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2 90° Corner Connector	82
C2399. Pass-Through Cover	75
C2410. Rail-Hanging Work Surface	74
C2512. Sink Module	106
C2513. Sink Unit for Storage	107
C2520. ADA Sink Support, Systems Mounted	108
C2910. Work Surface Backsplash	78
C2920. Backsplash Filler	79
C3011. B-Style Shelf	135
C3019. Add-On Shelf with Lip	137
C3091. B-Style Shelf with Sloped Top	136
C4310. Workstation, Electric Height Adjustment	85
C4311. Workstation, Manual Height Adjustment	87
C4312. Work Table, Electric Height Adjustment	89
C4313. Process Table, Electric Height Adjustment	83
C4320. Riser Frame Divider	101
C4321. Riser Frame End Cover	102
C4322. Power Strip	103
C4323. Marker Board	104
C4324. Flat Panel Mount	105
C4326. Reagent Shelf	84
C4327.	
C4700. Work Table Surface, Rectangular	97
C4900. Work Table Reagent Shelf	100
C5009. Interior Filler	112
C5010. Storage Unit	109
C5014. Storage Unit with Sliding Doors	121
C5015. Storage Unit Shelf	110
C5016. Storage Unit Door	111
C5017.	
C5020. Shelf Unit	122
C5021.	
C5022. Shelf Unit, Diagonal Corner	124
C5024. Shelf Unit, 90° Corner	125

# Index by Product Number *continued*

Index: Product Number

C5025.	Add-On Shelf, Shelf Unit	126
C5026.	Shelf Unit Back Panel	128
C5027.	Shelf Unit Door	130
C5028.		
C5032.	Shelf Unit Double Doors	131
C5036.	Shelf Unit Wall Attachment Panel	133
C5040.	Storage Case	114
C5041.	Storage Case Shelf	115
C5042.	Storage Case Pullout Shelf	116
C5043.	Storage Case Drawer	117
C5044.		
C5045.	Storage Case Bottle Drawer	118
C5048.	Lockable Drawer Bearer	113
C5050.	Storage Case Leg	120
C5051.	Light Bracket	129
C5052.	Storage Case Door	119
C5053.		
CO201	C Lid	160
CO203	C Tray/Shelf	161
CO204	CS Tray/Shelf	157
CO205	C Frame	158
CO207	Drawer	165
CO208		
CO209		
CO212	L Cart	195
CO215.	Draw Rod, Seismic	12
CO219	Locker Tambour Door	150
CO228	Drawer Bearer Package	175
CO231	Subcontainer	181
CO232		
CO233		
CO234		
CO235		
CO236		
CO238		
CO240		
CO242	Subdivider Vane	162
CO252	C Wire Shelf	163
CO266	Gripper Package	164
CO274	Label Clip	177
CO277	Locker Lock Kit	153
CO279	Locker Lock Kit, Keyless	155
CO282	Standard Rail	5
CO284	Rail Hardware	6
CO285		
CO290	Locker Hinged Door	151
CO338	End Cap for Standard or Seismic Rail	10

CO365	Dispensing Rail	179
CO366		
CO367		
CO368		
CO369	Subcontainer	181
CO370	Drawer with Lockable Lid	168
CO371	Locked Drawer with Lockable Lid	170
CO372	Drawer with Lockable Lid	168
CO373	Locked Drawer with Lockable Lid	170
CO374	Drawer with Lockable Lid	168
CO375	Locked Drawer with Lockable Lid	170
CO376	Subcontainer Handle	186
CO377		
CO378	Subcontainer Divider	185
CO379		
CO380		
CO381	Subcontainer Label Flag	188
CO382		
CO396	Lockable Drawer Bearer	176
CO397.	Sloped Shelf Unit Top	132
CO398	Shelf Label Clip	127, 138
CO459	Subcontainer Handle Label	187
CO462	Label Clip	177
CO469.	Pullout Catheter Rack	189
CO479	Wall Rail, Seismic	11
CO481	Floor Anchor Bracket, Action Office® Series 1	14
CO485	Floor Anchor Bracket, Action Office® Series 2	15
CO487.	Wall Strip Cover	19
CO508.	Clear Lid	174
CO521.	B-Style Rail-Hanging Shelf	134
CO540	A-Size Drawer Organizer Kit	183
CO541		
CO542		
CO543		
CO544		
CO545		
CO546	B-Size Drawer Organizer Kit	184
CO547		
CO552	Floor Anchor Adapter Kit	17
CO553	Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	18
CO555	Counter Top	159
CO556		
CO561	Locker, with Tambour Door	139
CO562	Locker	143
CO563	Locker on Wheel Base	149
CO564	Locker on Wheel Base, with Tambour Door	144
CO565	Cassette Assembly	178



CO568	Floor Anchor Bracket, Co/Struc® Lab Module	13
CO570	Locker, with Hinged Door	141
CO574	Locker on Wheel Base, with Hinged Door	147
CO600	Laptop Cart	190
CO601	Flat Panel Cart	191
CO602	Coiled Power Cord	192
CO605	Utility Basket	193
CO608	Cable Manager	194





## 20-Day or Less/Assigned Lead-Time Order Information — Healthcare

GSA

### 20-Day or Less/Assigned Lead-Time Programs

20-day or less and assigned lead-time products are designed for established Herman Miller dealers with pre-approved credit.

Products are distributed from multiple manufacturing facilities across the United States. All products and options not designated by an Assigned Lead-Time icon  will ship in 20 business days or less after being acknowledged by Herman Miller. Products and options designated by an  are on the assigned lead-time program and will ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

### Ordering Procedure

Please place orders through Order Manager. If this is not available to you, place orders through mail or fax. Verbal purchase orders will not be accepted.

Order Entry fax number: (616) 654 3085.

For more information, contact your Customer Care representative at: (866) 854 3048 ext 3400

### Shipments and Delivery

Per Herman Miller Terms and Conditions.

### Changes and Cancellation

Per Herman Miller Terms and Conditions.



# Packaging/Keyed-Alike Information — Co/Struc® System

## GSA

### Packaging

Packaging varies according to individual product. A surcharge may be added for special packaging requests. Contact Customer Care for specific information.

### Keyed-Alike Locks

All components with locks may be specified keyed alike.

Keyed-alike Co/Struc components with a “CO” or “CV” prefix will be shipped with lock plugs installed. Mobile Technology Support products (CO600, CO601, and CO606) cannot be ordered keyed alike.

Keyed-alike Co/Struc components with a “C” prefix will be shipped with lock chassis installed. Lock plugs will be shipped separately, tagged for product and location as specified. Installation of lock plugs will be on-site according to instructions provided. Herman Miller, Inc., will not be responsible for lock plug installation.

To order keyed-alike locks for “C” prefix components, specify keyed alike. Then place a separate order for plugs and locks:

1. List plug part number and key number (232092-XXX); select a key number (or numbers) between 226 and 324.
2. List quantity of locks per key number.
3. Indicate this is a no charge order.

Keyed-alike locks are field installed.

A Medeco master key can be used to open any lock on Co/Struc components with a “CO” or “CV” prefix, except Mobile Technology Support products CO600, CO601, and CO606. Customers interested in ordering master keys should contact Customer Care with a written request. To remove lock plugs from existing components, loosen the hex nut inside the product and remove the plug. These locks are available in key numbers 01 through 99. To remove lock plugs from locker with electronic lock, remove the cover from the e-lock mounting plate (be careful not to separate battery wires from unit). Next, loosen the hex nut inside the cover and remove the plug.

A Hudson master key can be used to open any lock on Co/Struc components with a “C” prefix. To remove lock plugs from existing components, a removal key or lock change tool is required. Customers interested in ordering master keys, removal keys, or lock change tools should contact Customer Care with a written request. Locks for “C” prefix components are available in key numbers 226 through 324 but are ordered as key numbers 01 through 99. Use the following chart to cross-reference the number selection.

“C” lock order numbers	Actual “C” lock plug numbers
1	226
2	227
3	228
4	229
5	230
6	231
7	232
8	233
9	234
10	235
11	236
12	237
13	238
14	239
15	240
16	241
17	242
18	243
19	244
20	245
21	246
22	247
23	248
24	249
25	250
26	251
27	252
28	253
29	254
30	255
31	256
32	257
33	258
34	259
35	260
36	261
37	262
38	263
39	264
40	265
41	266
42	267
43	268
44	269
45	270
46	271
47	272
48	273
49	274
50	275

“C” lock order numbers	Actual “C” lock plug numbers
51	276
52	277
53	278
54	279
55	280
56	281
57	282
58	283
59	284
60	285
61	286
62	287
63	288
64	289
65	290
66	291
67	292
68	293
69	294
70	295
71	296
72	297
73	298
74	299
75	300
76	301
77	302
78	303
79	304
80	305
81	306
82	307
83	308
84	309
85	310
86	311
87	312
88	313
89	314
90	315
91	316
92	317
93	318
94	319
95	320
96	321
97	322
98	323
99	324

For more information about our products and services or to see a list of dealers, please visit us at [HermanMiller.com](http://HermanMiller.com) or call (800) 851 1196.

© 2025 Herman Miller, Inc., Zeeland, Michigan


#### Distribution Rights

Working together since 1956, Herman Miller and Vitra are the authentic sources of designs by Charles and Ray Eames, Alexander Girard, and George Nelson. All the designs shown in this book are protected by law. Herman Miller has obtained all rights to make and sell these designs. Also, Herman Miller holds exclusive worldwide manufacturing and distribution rights for these products with the following exceptions:

Charles and Ray Eames, George Nelson, and Alexander Girard—Worldwide distribution rights for furniture in all areas except Europe and the Middle East. For those areas, please contact Vitra.

Isamu Noguchi—Worldwide distribution rights for the Noguchi Rudder Table. Distribution rights for the Noguchi Table in North America and Central America only. For all other areas, please contact Vitra.

Magis—Exclusive distribution rights in the US and Canada. For all other areas, please contact Magis.

® HermanMiller, , Action Office, Aeron, Aside, Canvas Office Landscape, Caper, Cosm, Co/Struc, Eames, Embody, Ethospace, Exclave, Fine-Tune, FlexFront, Formcoat, Ground Cloth, Intent, Layout Studio, Lino, Maharam, Meridian, Mirra, Motia, Nelson, Pellicle, Perspectives, Plex, PostureFit, PostureFit SL, Sayl, Setu, Tu, Verus, Y-Tower and Zeph are among the registered trademarks of Herman Miller, Inc., and its owned subsidiaries.

™ AireWeave, Bubbletack, Canvas Vista, Cellular Suspension, Commend, Compass, Connect, DOT, Durawrap, Everywhere, Flex-Edge, Headway, Hopsak 2, Iota, Keyless Entry, Latitude, Lyris 2, Mora, Multiscrim, Nevi, OE1, OE1 Workspace Collection, Overlay, Renew, Stackable, Swoop, Thrive, TriFlex, Twist and Valor are among the trademarks of Herman Miller, Inc., and its owned subsidiaries.

Chemsurf® is a registered trademark of Wilsonart International.

Corian® is a registered trademark of E.I. du Pont de Nemours and Company Corporation.

FLEXNET™ is a trademark of Milliken & Company.

Floor Saver Glide® is a registered trademark of Carpin Manufacturing, Inc.

Forest Stewardship Council® is a registered trademark of the Forest Stewardship Council A.C. (FSC). We are FSC® certified (FSC® C102895).

FSC® is a registered trademark of the Forest Stewardship Council A.C. (FSC).

MicrobeCare™ is a trademark of Parasol Medical LLC.

Microsoft® and Natural® Keyboard are registered trademarks of Microsoft Corporation.

Place Textiles® is a registered trademark of Place Textiles.