



**The beautiful** design of Restore extends from the waiting area to patient rooms with the patient recliner and sleeper sofa.

**The customizable** recliner is available with a variety of options, including a folding tablet, casters and push bar, and IV pole.

**Three** seating positions are featured; the full recline position may be specified with the trendelenberg mechanism.

**Convert** the sleeper sofa into a bed with minimal effort by folding down the back cushion, which is secured to ensure separate seating and sleeping surfaces.

**Stow** items in the storage compartment by lifting the seat, which incorporates a lid stay hinge a weightless feel and to allow the lid to remain where positioned.

