

# Chrys

Especially in today's world, creating products and using materials that can be properly cleaned to protect the health and wellbeing of others is of paramount importance. This is why BIFMA created strict cleanability guidelines for the entire industry. Here at OFS and Carolina, however, we decided to take our own personal cleaning guidelines even further. We commit to crafting products that can withstand stronger cleaning agents for a deeper clean and for better-protected environments for the people you love.

	Ammonia (16:1)	Bleach (full strength)	Hydrogen peroxide	3M TB quat - straight	Purell hand sanitizer	Cavicide	Iso alcohol	BAC-50	Sani-cloth HB	TSP	Wex-cide RTU	Virox RTU	Betadine	Oxivir TB	Betco pH7Q	Virex TB	Oxycide
Agave	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs
Jute	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs
Quill	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs
Sea Salt	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs	24 hrs

\*BIFMA's guidelines provide structure for furniture manufacturers and healthcare professionals to understand typical cleaners, disinfectants, cleaning methods, and expectations for recommended furniture performance/standards when exposed to these cleaners. To test performance, BIFMA requires a one-inch square of the solution to be applied to the surface for 15 consecutive minutes. After, it can be cleaned with distilled water and dry for 1 hour before it is tested for surface discoloration. If there is discoloration, it fails the test.